

ECLIPSE

Suspension Sleeve

The NEW Eclipse Suspension Sleeve provides extreme durability and comfort for your patient. Imagine not having to replace suspension sleeves every 3 months. A unique combination of materials and smart design, ST&G has created a superior suspension sleeve that offers high performance and secure fit.



- Reinforced silicone and strategically located Parylene® coating provides maximum durability.
- Dual sleeve design with non-slip inner strip of surface provides a maximum airtight seal which eliminates slippage.
- This revolutionary patent pending design helps keep the sleeve in place even during the most strenuous activity.
- Depending on activity, the outer Eclipse sleeve easily folds down over itself to reduce the length.
- The Parylene® coating does not inhibit knee flexion which allows for excellent performance and comfort.
- The Eclipse offers a platinum cured medical grade silicone with a Parylene® surface which reduces surface friction between the sleeve and skin helping in the prevention of skin shear for fragile skin and blister reduction.
- The smooth and silky coating allows for clothing to naturally hang and not bunch up.
- Warranty: 120 day warranty against defect and workmanship.

P/N	Description	Size	Price
3400B-S	BLACK Suspension Sleeve / Small	Small (41): 41cm – 49cm	\$197.00
3400B-M	BLACK Suspension Sleeve / Medium	Medium (49): 49cm – 55cm	\$197.00
3400B-L	BLACK Suspension Sleeve / Large	Large (55): 55cm – 62cm	\$197.00





ECLIPSE

Suspension Sleeve Instructions



The Eclipse is a unique two seal suspension system for the active amputee. The Eclipse is designed to be used in variable setups.



The outer Eclipse sleeve can be rolled or folded down, using only the inner sleeve to suspend. For best suspension, use both sleeves in conjunction with each other. Many variations are being employed by the amputees currently wearing the Eclipse.



Shiny strip of silicone

Donning the Eclipse



1. Don the short inner sleeve on a clean socket low enough that the shiny end is completely sealed against the socket wall. Be sure the proximal end of the inner sleeve covers the medial and lateral ears.

Shiny strip of silicone



2. Don the outer sleeve so that the edge of the outer sleeve is in line with the edge of the inner sleeve. Make sure the shiny strip of silicone of the outer sleeve is touching the inner sleeve's shiny silicone. Once in position roll the outer and inner sleeve down.



3. Don the prosthesis.



4. Pull up the inner sleeve.



5. Roll up the outer sleeve. Use extra care to force out any trapped air.



Durability suggestions:

1. Trim the socket edges so there is no extra reverse flairs.
2. Carbon sockets need a coating (super glue) applied to all edges. Carbon is a strong abrasive.
3. Roll the prosthetic socks over the ears before donning the inner Eclipse sleeve.
4. When a hole occurs rotate the outer or inner Eclipse sleeve 5 degrees. The rotation will seal the leak.

