

	PRODUCT	MOBILITY SCALE	Articulating Ankle 	Heel Height Adjustable 	Vertical Shock 	Multi-Axial 	Split Toe 	Sandal Toe 	SIZES	BUILD HEIGHT	MAXIMUM BODY WEIGHT	PRODUCT WEIGHT(AVG)	FREEDOM PROSTHETIC FEET
	Renegade® (RS1)				✓			✓	22-31cm	6 3/8" - 7 7/8"	365 lbs.	765 g	
	Renegade® AT (R11)				✓	✓	✓		22-31cm	6 3/8" - 7 7/8"	365 lbs.	765 g	
	Renegade® LP (RS6)				✓			✓	22-31cm	5" - 5 5/8"	365 lbs.	765 g	
	Renegade LP-AT (R16)				✓	✓	✓		22-31cm	5" - 5 5/8"	365 lbs.	765 g	
	Thrive® (FS5)						✓	✓	22-31cm	5 3/4" - 7"	365 lbs.	765 g	
	Highlander® (FS3) / MAX (FS3-H5)						✓	✓	22-31cm	4 1/4" - 5 1/2"	500 lbs.	595 g	
	Kinterra® (ROM)		✓			✓	✓	✓	22-31cm	4 3/4" - 4 7/8"	275 lbs.	652 g	
	Agilix				✓	✓	✓	✓	22-31cm	6 1/8" - 7 1/2"	365 lbs.	635 g	
	Pacifica® LP (FS4)						✓	✓	22-31cm	1 7/8" - 2 3/8"	365 lbs.	570 g	
	Promenade® (VS2)						✓		19-31cm	4 5/8" - 5 3/8"	325 lbs.	595 g	
	Senator® (VS1)								22-31cm	5 1/8" - 6"	300 lbs.	652 g	
	Sierra® (FS1)					✓	✓	✓	22-31cm	5 3/4" - 7 3/8"	365 lbs.	652 g	
	DynAdapt® (F10)					✓	✓	✓	22-31cm	6 1/8" - 7 1/8"	365 lbs.	623 g	
	Runway® (RS2) / HX (RS2-K1)			✓				✓	22-28cm	3 3/8" - 3 7/8"	255 lbs.	26 oz.	
	Silhouette® (RS4/8)				✓	✓			22-31cm	7" - 7 7/8"	365 lbs.	652 g	
	Silhouette® LP (RS5/9)				✓	✓			22-31cm	6 5/8" - 7 1/2"	365 lbs.	538 g	
	Silhouette® LP-VS (R15/R19)				✓	✓			22-31cm	6 5/8" - 7 1/2"	365 lbs.	538 g	
	Silhouette® VS (R14/R18)				✓	✓			22-31cm	7" - 7 7/8"	365 lbs.	652 g	
	Defender™ Pediatric (FS3)								16-21cm	2 3/8"-3"	100 lbs.	226 g	
	LP Symes™ (LP2)								22-31cm	1 3/4" - 2 1/8"	365 lbs.	538 g	
	Catapult® Run (FX6) / Sprint (FX7)								N/A	9 1/2"	365 lbs.	510 g	
	Freestyle™ Swim (LP2-W2)								22-31cm	4" - 4 3/8"	255 lbs.	793 g	
	Nitro™ LP Run (FS9)								N/A	7 1/4"	365 lbs.	311 g	
	Nitro™ Run (FS6)								N/A	9 1/2"	365 lbs.	283 g	
	Slalom™ Ski (SKI)								custom	5 1/2"	250 lbs.	1050 g	
	WalkTek® (VS3)					✓		✓	21-30cm	3 3/8" - 3 3/4"	300 lbs.	623 g	

K2

3 Skin Tones Available

