

# Custom Order Form: Knee Supports

Please check all boxes that apply

## Material

- Neoprene
- Coolprene (perforated Neoprene)
- Coolcel (breathable, latex-free)

## Thickness and Colour

(Neoprene + Coolprene only)

- |                              |                                |
|------------------------------|--------------------------------|
| <input type="checkbox"/> 2mm | <input type="checkbox"/> Black |
| <input type="checkbox"/> 3mm | <input type="checkbox"/> Grey  |
| <input type="checkbox"/> 6mm | <input type="checkbox"/> Blue  |

## Hinges / Stabilizers

- Polycentric hinges (32 style)
- Multistop hinges (39 style)
- Spiral stays (34 Style)  1  2  3

## Straps

- Elastic circumferential (32 Style -1 superior, 1 inferior)
- Velcro medial to lateral (34 Style - 1 superior, 1 inferior)
- Velcro circumferential (39 style - 2 superior, 2 inferior)


## Additional Options

- Wrap around closure (no buttress)
- Open top closure
- Popliteal diamond cut out
- Pull Tabs
- Other: \_\_\_\_\_

## Company Information

Company Name : \_\_\_\_\_ Account #: \_\_\_\_\_  
 Patient Name : \_\_\_\_\_  
 Contact Name : \_\_\_\_\_ Phone # \_\_\_\_\_  
 Order Date : \_\_\_\_\_ P.O. # \_\_\_\_\_

Please mark all measurements in inches. Please supply two additional measurement half way up and down the brace if total brace length exceeds 12".



Distance from knee center to top of brace (normally 6")

Distance from knee center to bottom of brace (normally 6")

Above Patella Circumference

Knee Center Circumference

Below Patella Circumference

Total Length

Left Leg

Right Leg

## Buttress Type

- Foam  Removable
- Sewn in superior
- Sewn in inferior
- Tube  Donut
- Lateral J
- Split J
- None (patella hole only)