

Ortho Active is now distributing products from Pavis Italy. Pavis manufactures high quality orthopedic supports with some unique designs and premium materials. The Pavis New Edge support products are all made with a 100% velour cotton lining for comfortable skin contact and a Velcro-sensitive exterior for easy application.

6622 New Edge Hinge Knee Wrap

100% Breathable

A silicone ring provides patella stability, and removable, mouldable polycentric hinges provide medial and lateral support to the MCL ligaments of the knee. The offset contoured hinge should be inserted on the medial side while the straight hinge is for the lateral side. A lycra popliteal area provides good circulation. This brace opens completely for easy application and has two elastic straps to adjust the compression and increase support.

Size	Thigh Circ. (10 cm above knee centre)
Small	35-40 cm
Regular	40-48 cm
Extra	48-58 cm



6633 New Edge Wrist Lacer

100% Breathable

The three Velcro closures with adjustable anchor point for the laces and ambidextrous design make this a versatile wrist brace. Three mouldable removable stays provide excellent support for the wrist. The contour stay inserted on the palm side and the straight stay inserted on the dorsal side determine if it is a left or right brace. The brace opens nicely, making it easy to apply, and the lace/Velcro closure system is easy to adjust and tighten.

Size	Wrist Circumference
Small	13.5-15.5 cm
Regular	15.5-18.5 cm
Extra	18.5-21 cm

6637 New Edge Wrist Thumb Lacer

100% Breathable

The three Velcro closures with adjustable anchor point for the laces make this an easy-to-fit wrist/thumb brace. Four mouldable, removable, pre-contoured stays provide excellent support for the wrist and thumb and are very comfortable. The brace opens nicely, making it easy to apply, and the lace/Velcro closure system is easy to adjust and tighten. Specify Left or Right.

Size	Wrist Circumference
Small	10-14 cm
Regular	14-20 cm
Extra	20-26 cm



6639 New Edge Elbow Lacer

100% Breathable

A double wave shaped pressure pad provides effective pain reduction. A moveable anchor point for the lace closure system makes it easy to fit and adjust the size of this brace. Simply place the pad of the strap on the forearm muscle, apply the Velcro anchor to the outside of the support, and tighten the lace/Velcro strap for easy application and tension adjustment.

Size	Forearm Circ.
Regular	22-32 cm
Extra	30-45 cm

6640 New Edge Patella Stabilizer

100% Breathable

Two removable patella C-shaped buttresses that work in combination with the central straps allow optimization of the position of the buttresses and support for the patella. Removable, lightweight, flexible lateral stays improve knee stability. Lycra popliteal area allows for maximum circulation. Lace elastic strap closure system provides optimum adjustability and comfort.

Size	Thigh Circ (10 cm above knee centre)
S	32-36 cm
M	36-40 cm
L	40-45 cm
XL	44-49 cm





6655 Wellness Boxer Hernia Support

Bacteriostatic • 100% Breathable

The Boxer Hernia Support integrates the design of a traditional hernia support with high-strength elastic and removable hernia support pads into a stylish and comfortable pair of boxer shorts. The hernia pads are anatomically shaped to reduce the hernia for pain relief. High-tech athletic materials are used to provide breathability, antibacterial action and high durability for full-time wear and washing. Nanotechnology yarn ensures the bacteriostatic quality remains unchanged for more than fifty washing cycles. The boxer style is very effective for active use during sport activities.

Size	Pelvis Circ.
S	71-82 cm
M	82-94 cm
L	94-105 cm
XL	105-117 cm
XXL	117-130 cm

6656 Wellness Brief Hernia Support

Bacteriostatic • 100% Breathable

This brief-style hernia support integrates the design of a traditional hernia support with high strength elastic and removable hernia support pads into a stylish and comfortable pair of briefs. The hernia pads are anatomically shaped to reduce the hernia for pain relief. High-tech athletic materials are used to provide breathability, antibacterial action and high durability for full time wear and washing. Nanotechnology yarn ensures the bacteriostatic quality remains unchanged for more than fifty washing cycles. The boxer style is very effective for active use during sport activities.

Size	Pelvis Circ.
S	71-82 cm
M	82-94 cm
L	94-105 cm
XL	105-117 cm
XXL	117-130 cm



6678 Wellness Ostomy Support Belt

Bacteriostatic • 100% Breathable

A 9" wide ostomy support belt with adjustable elastic compression to prevent abdominal hernias. The colostomy hole is 7cm in diameter and can be adjusted up to 11 cm in diameter by cutting on the circular welding of the front panel material. 100% soft cotton fabric is in contact with the skin combined with Pavis Airflex woven technology including bacteriostatic and hypoallergenic fabric for maximum comfort and support. A removable bag cover in ultralight elastic cotton attaches to the front of the support.

Size	Waist Circ.
S	75-90 cm
M	90-115 cm
L	115-135 cm