



EASY WALK® Fixation

Easy Walk® Fixation is created to maximize ankle support and optimize biomechanical gait cycle. Strong lateral fixation of the lower leg provides maximal medial-lateral control of the foot-and-ankle complex. Due to rigidity of the support strut, the major energy return during the gait is carried out by the footplate.

INDICATIONS

- severe foot-drop*
- severe ankle instability*
- spasticity*
- foot and calf extensors after surgery*
- proprioceptive deficit*
- ankle arthritis*



Item No. Left	Item No. Right	Size	Shoe size	Height, cm
AFO10F=LXS	AFO10F=LXS	XS	34-36	32.5
AFO10F=LS	AFO10F=RS	S	36-39	33.5
AFO10F=LM	AFO10F=RM	M	39-42	34.5
AFO10F=LL	AFO10F=RL	L	42-45	35.5
AFO10F=LXL	AFO10F=RXL	XL	45-47	37.5