



Kite
AFO 

www.orthobaltic.eu
PATENT PENDING



Easy Walk® KITE

Ortho Baltic is proud to introduce the latest addition to the Easy Walk® range of carbon fibre AFOs, this time dedicated to the youngest patients – Easy Walk® KITE*.

Orthosis is intended for children with dorsiflexion weakness and is specifically designed to control medial-lateral movement and stabilize the ankle. At the same time, the carbon footplate assists the gait by flexing and providing the energy return for toe-off.

Durability of Easy Walk® KITE was proven by the tests, carried out by notified body in Berlin, in compliance with ISO 10328 standard. The orthosis successfully withstood over one million load cycles, as well as an ultimate strength test with max force of 3220N.



INDICATIONS

- Unstable or low tone gait
- Foot-drop
- Cerebral palsy
- Spina bifida



Item No. Left	Item No. Right	Size	Shoe size	Height, cm
AFO10K=LS	AFO10K=RS	S	20-24	21
AFO10K=LM	AFO10K=RM	M	24-27	25
AFO10K=LL	AFO10K=RL	L	27-30	29
AFO10K=LXL	AFO10K=RXL	XL	30-33	34

*Easy Walk Kite® Patent no.: LT6591