

## foot size & category chart

Applies to Forever 6", 7" and 9".

Body Weight		Sizes		
Weight in kg	Weight in lbs	23 - 24	25 - 26	27 - 29
35 - 50	77 - 110	Category 1	Category 1	N/A
50 - 65	110 - 143	Category 1.5	Category 1.5	Category 1.5
65 - 80	143 - 176	Category 2	Category 2	Category 2
80 - 95	176 - 210	Category 2.5	Category 2.5	Category 2.5
95 - 110	210 - 242	N/A	Category 3	Category 3
110 - 125	242 - 275	N/A	Category 3.5	Category 3.5

The Levitate Forever comes in size 23-29 and covers weight up to 125 kg / 275 lbs.

Depending on the users activity level you can go one category higher or lower.



## foot size & category chart

Applies to Forever LP

Body Weight		Sizes		
Weight in kg	Weight in lbs	23 - 24	25 - 26	27 - 29
35 - 50	77 - 110	Category 1	Category 1	N/A
50 - 65	110 - 143	Category 1.5	Category 1.5	Category 1.5
65 - 80	143 - 176	Category 2	Category 2	Category 2
80 - 95	176 - 210	Category 2.5	Category 2.5	Category 2.5
95 - 110	210 - 242	N/A	Category 3	Category 3
110 - 125	242 - 275	N/A	Category 3.5	Category 3.5

The Levitate Forever LP comes in size 23-29 and covers weight up to 125 kg / 275 lbs.

Depending on the users activity level you can go one category higher or lower.



foot size & category chart



The Levitate Blade Sereies comes in 2 different sizes and covers a wide range of user needs and covers weight up to 130 kg / 286 lbs.

Blade 8"

Body weight		Activity level	
Kg	Lbs	Low	High
10-15	0-33	Category 1	Category 2
15-30	33-66	Category 2	Category 3
30-45	66-99	Category 3	Category 4
45-60	99-132	Category 4	Category 5
60-75	132-165	Category 5	Category 6
75-90	165-198	Category 6	N/A

Blade 10"

Body weight		Activity level	
Kg	Lbs	Low	High
30 - 40	66 - 88	Category 2	Category 2.5
40 - 50	88-110	Category 2.5	Category 3
50 - 60	110 - 132	Category 3	Category 3.5
60 - 70	132 - 154	Category 3.5	Category 4
70 - 80	154 - 176	Category 4	Category 4.5
80 - 90	176 - 198	Category 4.5	Category 5
90 - 100	198 - 220	Category 5	Category 5.5
100 - 110	220 - 242	Category 5.5	Category 6
110 - 130	242 - 286	Category 6	N/A

## traction footshell

overview

Engineered for performance and everyday reliability, the Traction footshell features a streamlined, technical design that delivers enhanced grip, long-lasting durability, and waterproof protection—ensuring dependable function across a wide range of daily environments.

key features

**Exceptional grip:** Designed to be worn inside a shoe, Traction provides secure grip and responsive support during high-activity use.

**Technical aesthetic:** Made in a bold black finish and minimal design, Traction breaks away from the traditional look of footshells.

**Waterproof Design:** Built for wet environments with integrated drain holes that let water flow through easily—ideal for barefoot use at home, by the pool, on the terrace, or at the beach.

compatibility

Forever Series, Forever LP.



## human footshell

overview

The Human footshell is a robust and durable everyday option, designed to be worn comfortably inside a shoe. With a natural shape and reliable fit, it protects the prosthetic foot while providing smooth function for daily activities.

key features

**Low ankle cut:** Allows the footshell to fit seamlessly into a wide range of shoe styles, offering flexibility and comfort without compromising stability.

**Hardwearing material:** Built from durable, abrasion-resistant materials that extend the lifespan of the prosthetic by protecting it from daily wear, impacts, and environmental stress.

compatibility

Forever Series, Forever LP.

Human footshell - Black



Human footshell - Beige

