



The lifespan of Formthotics varies based on usage but typically ranges from 12 to 24 months. Factors such as the intensity and frequency of use, activity level, and individual foot characteristics can influence their durability. Proper care and maintenance also play a crucial role in extending their lifespan.

They are made from the unique, proprietary PE foam called Formax, a type of chemically cross-linked, closed-cell polyethylene foam. This foam is very different from EVA that most other foam insoles are made from. Formax is completely thermoformable (heat mouldable), which means they can change shape at relatively low temperatures. Its closed-cell structure means it is completely waterproof.

Formthotics products come with a 30-day satisfaction guarantee. Check with your distributor or retailer for specific warranty information.