

## Home Healthcare

Section ..... **A**



## Cervical

Section ..... **B**




## Elbow & Shoulder

Section ..... **C**



## Torso

Section ..... **D**



## Hip & Groin

Section ..... **E**



## Wrist & Hand

Section ..... **F**



## Leg & Knee

Section ..... **G**



## Ankle

Section ..... **H**



## Footcare

Section ..... **I**




## Pediatric

Section ..... **J**



## Order Forms

Section ..... **K**



### Contact Us

Phone: 604-520-3414

Toll Free Phone: 1-800-663-1254

Fax: 604-520-1193


Toll Free Fax: 1-888-558-6565


Email: [orders@orthoactive.com](mailto:orders@orthoactive.com)

For over 40 years, Ortho Active has strived to deliver the highest quality products and raise the standard for customer service with unwavering integrity. We thank you, our customers, who have supported us on the front lines of patient care by using our products to change people's lives for the better. We appreciate your business and consider it a privilege to serve you both now and for many years to come.

## Legend

Indicates product # if no order chart is present      Product name      Indicates manufacturer (if applicable)



**201100 Stabilo Hip Protector** 

This high quality hip protector is designed to protect the geriatric hip from fracture due to falling. Made from cotton and Lycra with a highly elastic waistband and leg borders. Bilateral protectors are made of shock absorbing visco-elastic soft foam and are easily removed. Select best fit size according to chart below.

Product #	Description
<b>201100</b>	Stabilo Hip Protector
<b>201120</b>	Replacement Pads
<b>201200</b>	Briefs Only
<b>201250</b>	Stabilo Hip Incontinence

**SIZING:**

Size	Waist Circ.	Hip Circ.
S	27" - 29"	35"
M	27" - 31"	37"
L	31" - 32"	38"
XL	32" - 34"	40"
XXL	34"+	42"+

Product image      Order chart with product descriptions      Product description      Sizing chart with required measurements

### Order Notes For Softgoods

**Retail Packaging** - Please place an "R" before the product number to receive retail packaging (i.e. R3193 instead of 3193). If another manufacturer is indicated, the item automatically comes in retail packaging if available.

Knee brace type (hinged)      Indicates Neoprene (non-perforated, standard material)

**Neoprene Products -**      i.e. **32 P/3**      Neoprene thickness (for knee braces only)

Buttress type (Universal felt)      Colour is black Neoprene with black plush lining

**Coolcel** products are manufactured with breathable, open cell foam and include only latex free materials. All Coolcel product numbers begin with 31; for example a Coolcel Hinged Knee Brace is item number 3132.

**Airflex** products are manufactured with moisture wicking, stretchy, and smooth material that is also latex free. All Airflex product numbers begin with 37; for example an Airflex Thumb Support is item number 3767.

## Fabric Types

Throughout our Ortho Active manufactured products, there is a range of fabric types to choose what will work best for your client's needs. Please find details below on the different fabrics and where they are best used.

### Neoprene

Neoprene is an extremely durable material, and is the standard fabric (if a selection is not otherwise made) for our Ortho Active products. Neoprene is a non-perforated latex-free foam material which provides the most warmth of the fabric alternatives. Ideal for patients that can benefit from added insulation such as those with arthritis. Available in 3mm black.

### Coolcel

Coolcel is Ortho Active's breathable, latex-free Neoprene alternative. Constructed of a unique nylon/open cell foam combination, Coolcel is designed to give compression and support while allowing increased breathability and wicking action. Ideal for sport use or patients susceptible to skin irritation from heat rash. Available in black.

### Airflex

Airflex is extremely breathable material for cool comfort during active use. It is specifically designed for air circulation, heat transfer, and moisture wicking. It is stretchy and smooth for the highest level of comfort. Any activity where excess heat or perspiration could be felt under the brace is an ideal application for this fabric. Available in black.

## Proud Supplier of:



and more!

## **Pricing, Payment, and Pricing Compliance**

All prices are in Canadian dollars (CAD) and exclude applicable taxes and freight. Prices are subject to change or withdrawal without notice. Prices prevailing on date of shipment apply. Quotations are valid for thirty (30) days unless otherwise stated. Orders become binding only upon written confirmation by Ortho Active. Any Customer purchase-order terms that conflict with these Terms are expressly rejected unless agreed to in writing by Ortho Active.

Payment is required in full prior to shipment, unless credit terms have been approved in writing. Accepted payment methods: Electronic Funds Transfer (EFT) or e-transfer to the account designated by Ortho Active. Credit-card payments are subject to a 3% processing fee. Credit terms may be granted upon request and are subject to credit approval.

Bulk orders, booking orders, and other large-volume orders are subject to a 50% down payment at the time of purchase-order processing. Such orders are non-cancellable once confirmed.

Overdue balances on approved credit accounts are subject to 2% per month (24% per annum) interest. Ortho Active reserves the right to suspend or withhold shipments for accounts with outstanding balances.

Where a Manufacturer's Suggested Retail Price (MSRP) or recommended resale price is provided by Ortho Active, the Customer agrees to adhere to such pricing guidance when offering products for resale, unless otherwise authorized in writing by Ortho Active. Ortho Active reserves the right to modify MSRP from time to time. Updated pricing will be communicated through official price lists, electronic communications, or reseller notices. Failure to comply with MSRP requirements may result in suspension of supply, cancellation of pending orders, or termination of reseller authorization.

The Customer shall advertise products at prices materially below the published MSRP; engage in pricing practices that may devalue the brand, disrupt market fairness, or materially affect other authorized resellers; or utilize online channels, digital ads, or promotional mechanisms to circumvent MSRP guidance.

## **Delivery, Freight, Risk, and Title**

Unless otherwise stated, all shipments are F.O.B. Ortho Active warehouse (Coquitlam, BC).

## **Shipping**

The Customer pays freight on orders from \$50 to \$499.99 CAD (before tax). Free shipping applies to orders \$500 CAD and above, excluding expedited and/or large, bulky, or heavy items as designated by Ortho Active. We will select carriers and methods balancing speed, cost, and reliability. Customer shipping preferences or routing instructions must be provided at the time of order. Ortho Active will make reasonable efforts to accommodate such requests but will not be responsible for higher freight costs or delays when no routing instructions are given. Risk of loss transfers to the Customer upon transfer of goods to the carrier. Title passes to the Customer upon full payment. The Customer is responsible for maintaining appropriate transit insurance.

## **Damage or Loss**

All shipments are sold on an F.O.B. (free on board) Coquitlam basis. Ortho Active packages all products carefully, but once goods are transferred to the carrier, Ortho Active cannot accept responsibility for transit damage or loss. The Customer must inspect shipments immediately upon receipt. Any damage, loss, or missing product claim must be filed within 10 days from shipment date to Ortho Active.

Any claim not within this period may not be acknowledged or processed. Ortho Active will assist by providing documentation as required but cannot assume liability for freight-carrier damages or delays.

## **Short Shipments or Errors**

Claims relating to missing items, shipping errors, or quantity discrepancies must be reported to Ortho Active within five (5) business days of receipt of goods. Claims received after this period may not be accepted.

## **Product: Information, Use, Suitability, and Handling**

Product pricing, descriptions, specifications, and availability are subject to change without notice. Ortho Active makes reasonable efforts to ensure accuracy but does not guarantee that catalogue or website information is error-free.

Products which are medical devices supplied by Ortho Active are intended solely for use by qualified healthcare professionals. The Customer is solely responsible for determining the suitability of any product for its intended clinical use or patient application. Products sold by Ortho Active that are not classified as medical devices may be used in accordance with manufacturer instructions and applicable professional standards. No product is warranted as fit for a particular patient or purpose beyond the manufacturer's stated warranty.

The Customer must handle, store, and maintain all products in accordance with manufacturer instructions, professional standards, and Health Canada regulations. Ortho Active is not liable for damage or defect resulting from improper storage, handling, or modification.

## **Custom Products**

Ortho Active may supply custom-made or modified products to meet specific customer or patient requirements. Custom products include, but are not limited to, items manufactured to non-standard dimensions, specifications, or materials, or products altered at the customer's request. All custom-made or special-order products are final sale, non-returnable, and non-refundable, except where a confirmed manufacturing defect exists under the applicable warranty. Custom products are covered under their respective manufacturer or distributor warranty, or for Ortho Active-branded products, under the Ortho Active Product Warranty Statement. Customers must verify design details, measurements, and specifications before order placement. Production of custom items will proceed only after written confirmation or approval by the Customer.

## **Returns and RMA Procedure**

All returns require prior authorization through the Ortho Active Return Material Authorization (RMA) process. Returns must be requested within 30 days of invoice date. Products must be unopened, undamaged, and in original packaging. Medical Devices that have been fitted and/or worn are not returnable based upon Health Canada regulations on sterility. Returns are for credit on account only; no cash refunds. The Customer is responsible for return shipping costs and must provide a tracking number. Unauthorized or late returns may be refused. All returns are subject to a restocking fee of 20%.

For suspected product defects or device issues, the Customer must contact Ortho Active immediately before returning the product. Complaint related returns follow a separate Complaint Handling Process in compliance with the Medical Devices Regulations (SOR/98-282). Warranty returns will be processed per manufacturer terms.

## **Warranty**

Products are covered under their respective manufacturer or distributor warranties, as listed on the Manufacturer & Distributor Warranty Reference Guide, available at [www.orthoactive.com/warranty](http://www.orthoactive.com/warranty). Ortho Active is the first point of contact for warranty claims on all products purchased through Ortho Active. The warranty for Ortho Active-branded products is defined separately in the Ortho Active Product Warranty Statement.

Full Terms and Conditions can be found on our website: [www.orthoactive.com/terms-and-conditions](http://www.orthoactive.com/terms-and-conditions)

Product No.	Product	Page
A15 - L44	AquaShield	A2
80695/80696	CoolXChange® Compression & Cooling Gel Bandage	A2



## AquaShield

This simple, one piece, completely waterproof, orthopedic cast cover is reusable and offers a skid safe sole. Comes in different lengths/sizes. Strapless seal prevents water from entering and allows full submersion of arm or leg casts. For use with arm or leg plaster casts or to protect prostheses. Comes in pediatric sizes as well.

Product #	Width	Product #	Width
<b>A15</b>	Infant Arm 5.5" - 7.5"	<b>L15</b>	Infant Leg 7" - 9"
<b>A17</b>	Half Arm 17" (child)	<b>L20</b>	Half Leg 20" (child)
<b>A26</b>	Full Arm 26" (child)	<b>L33</b>	Full Leg 33" (child)
<b>A21</b>	Half Arm 21" (adult)	<b>L25</b>	Half Leg 25" (adult)
<b>A31</b>	Full Arm 31" (adult)	<b>L44</b>	Full Leg 44" (adult)



## CoolXChange®

R.I.C.E (Rest, Ice, Compression, Elevation) is a common first aid treatment for muscle strains and sprains, and may aid in quickening the recovery process. CoolXChange® advances R.I.C.E methodology by combining the cooling and compression stages, and is designed to provide a convenient and effective solution for pain relief. The CoolXChange® bandage assists in prevention, management and recovery from muscle pain, swelling and inflammation. It works instantly and is effective for up to 2 hours, decreasing skin temperature by over 10°C /18°F without causing ice burns.

Features:

- Ready to use anywhere, anytime
- Doesn't require refrigeration
- Provides immediate and continuous cooling for up to 2 hours
- Self-adhesive
- No mess, no smell
- Reusable



Product #	Description
<b>80695</b>	Regular - 3 x 47 in (ideal for ankles, feet, wrists, elbows and fingers)
<b>80696</b>	Large - 3 x 118 in (ideal for knees, hamstring, back, shoulder, or cut to fit)

Product No.	Product	Page
<b>Cervical Collars</b>		
<b>Foam Collars</b>		
3363	Cervical Collar w/ Stability Strap	B3
2020	Foam Cervical Collar (Beige)	B3
4506/4508/4510	Foam Cervical Collar (White)	B3

### 3363 Cervical Collar with Stability Strap (3" wide)



Contoured shape and vinyl stabilizer strap provide a comfortable and supportive cervical collar. 5/8" thick, medium firm foam is covered with a stockinette. The stabilizer strap is tan vinyl and closes with Velcro®. Foam length is 19" and stabilizer strap is 22" long.

**SIZING:** Universal (fits neck circumference 12" - 19")

### 2020 Foam Cervical Collar



Lightweight, medium density foam in a contoured design with both stockinette cover and tapered ends. Hook and loop closure. Provides warmth and mild mechanical restraint to the cervical region. Comes in beige.

**SIZING:**

Size	Neck Circumference	Brace Height
S	13" - 15"	3.5"
M	15" - 17"	4"
L	17" - 20"	4.5"

### Foam Cervical Collar



This cervical collar is made of breathable, high-quality foam padding to immobilize and support the neck as it heals. Soft, lightweight, and easy to wear. Has a slight contour to comfortably fit under the chin. Comes in 3 different heights for longer or shorter necks: 6 cm, 8 cm and 10 cm.

**SIZING:** Sizing: Measure neck height from chin to collarbone to choose the correct brace height, then measure neck circumference

#### 4506 (6cm)

S	25-30 cm
M	30-35 cm
L	35-40 cm
XL	40-45 cm
XXL	45-50 cm

#### 4508 (8cm) and 4510 (10cm)

S	35-40 cm
M	40-45 cm
L	45-50 cm
XL	50-55 cm
XXL	55-60 cm



Product No.	Product	Page
-------------	---------	------

## Elbow Supports

### Elbow Straps

61/62	Tennis Elbow Strap	C2
3162	Coolcel Tennis Elbow Strap	C2
3762	Airflex Tennis Elbow Strap	C2
6061	Tennis Elbow Strap	C2
022500	Bort EpiContur with Pad	C2
6639	Elbow Lacer	C3

### Elbow Sleeves

63	Tennis Elbow Sleeve	C3
3163	Coolcel Tennis Elbow Sleeve	C3
5564	Elastic Elbow Support	C3
5563	Dynamic Tennis Elbow Support	C3
055100	Bort Cubital Pad Elbow Sleeve	C4
122700	Bort Select EpiPlus Elbow Support	C4
816	Sportec Arm Compression Sleeve	C4

### Elbow Splints

5860	Ultra Cubital Tunnel Elbow	C5
------	----------------------------	----

Product No.	Product	Page
-------------	---------	------

## Arm Slings and Immobilizers

### Arm Sling/Shoulder Immobilizers

6011	Premium Shoulder Immobilizer	C4
531	SA Sling	C5
23191	Shoulder Abduction Pillow with Sling	C5

### Shoulder Suspension

540	Shoulder Suspension Brace	C5
-----	---------------------------	----



### 61/62 Tennis/Pickleball Elbow Strap

Tennis elbow strap with adjustable Velcro® closure to apply local pressure to the muscle of the forearm. Designed to provide pain relief by reducing muscle strain.

**Indications:**

- Tennis Elbow (Medial Epicondylitis)
- Golfer's Elbow (Lateral Epicondylitis)

**Fabric Alternatives:**

- Neoprene (available in black)
- Coolcel (with freezable pad for cold therapy)
- Airflex

Product #	Description
61	Neoprene without Tendon Pad
62	Neoprene with Tendon Pad
3162	Coolcel with Silicone pad
3762	Airflex without Tendon Pad

**SIZING:** See sizing chart at the bottom of the page



### 6061 Universal Tennis Elbow Strap

This support is ideal for sprains of the elbow ligament. The contact closure strap provides adjustable support and compression as you need it. The soft, latex-free Neoprene blend retains body heat for increased circulation to help relieve pain, promote healing, and keep your joint flexible. Wrap-around design makes it easy to don and doff.

**SIZING:** Universal



### 022500 Bort EpiContur with Pad

Epicondylitis soft brace for relieving the radial or ulnar muscle. A repositionable pad provides pressure, and the fitted cut enables a slip-free fit. Comes in a universal size, and is adjustable with a Velcro fastening. Designed to provide relief from chronic or recurrent epicondylopathy, tendopathies, overstraining, and tennis/golf elbow.

**SIZING:** Universal, fits forearm circumferences from 21-35cm

Product #	Description
022500	EpiContour brace
022500P	Silicone replacement pad

**SIZING:** Measure around largest part of forearm

Size	Forearm Circumference	Size	Forearm Circumference
XS	8" - 9"	L	11" - 12"
S	9" - 10"	XL	12" - 13"
M	10" - 11"	XXL	13" - 14"

### 6639 Elbow Lacer

A double wave shaped pressure pad provides effective pain reduction. A moveable anchor point for the lace closure system makes it easy to fit and adjust the size of this brace. Simply place the pad of the strap on the forearm muscle, apply the Velcro anchor to the outside of the support, and tighten the lace/Velcro strap for easy application and tension adjustment.

**SIZING:**

Size	Forearm Circ.
Regular	22-32 cm
Extra	30-45 cm



### Tennis Elbow Sleeve

Elbow sleeve with tennis elbow strap. To apply compression and warmth to the entire elbow area, as well as local pressure to the muscle of the forearm. Available with tendon pad.

Available in two fabric variations dependent on customer needs:

- **Neoprene:** unperforated, best for keeping in body heat
- **Coolcel:** Ortho Active's breathable, latex-free alternative ideal for sport use or patients susceptible to skin irritation from heat rash.

Product #	Description
<b>63</b>	Neoprene with strap
<b>63T</b>	Neoprene with tendon pad
<b>3163</b>	Coolcel with strap
<b>3163T</b>	Coolcel with tendon pad

**SIZING:** See chart on left page



### 5564 Elastic Elbow Support

**Features:**

- 3D knitting technology to provide superior fit and comfort
- Medical compression of 20-30mmHg
- Compression to improve circulation, reduce swelling, and support the active elbow for pain and tendon strains
- Knitted construction is breathable for cool comfort

**SIZING:**

Size	Elbow Circumference
S	8" - 9"
M	9" - 10"
L	10" - 11"
XL	11" - 12.5"
XXL	12.5" - 14"



### 5563 Dynamic Tennis Elbow Support

**Features:**

- Designed to provide warmth and compression to entire elbow
- Comfortable, breathable elastic will not overheat in warm weather or during activity
- Elastic strap supports the forearm muscle to promote healing of the tendonitis condition, and stretches to remain comfortable
- Two gel pads provide localized pressure to the forearm muscles, with dimpled gentle massaging surface
- Easily reverse the direction of elbow strap for application to either the right or left elbow

**SIZING:**

Size	Elbow Circ.
S	20 - 23 cm
M	23 - 26 cm
L	26 - 29 cm
XL	29 - 32 cm





## 055100 Bort Cubital Pad Elbow Sleeve

**Features:**

- Extensive padding and protection of entire elbow
- SOFTflex knitted fabric for more elasticity, stability and very comfortable wear
- Extensive silicone padding with recess at olecranon for dimensional stability during flexion
- Compression class II
- Seamless comfort zone in the antecubital fossa

**SIZING:** Forearm circumference

Size	Elbow Circ.
S	24 - 26 cm
M	26 - 28 cm
L	28 - 30 cm
XL	30 - 32 cm



## 122700 Bort Select EpiPlus Elbow Support



**Features:**

- Circular knitted elastic elbow for superior fit and comfort.
- Medial and lateral silicone pads provide intermittent massage and pressure
- The pressure on the pads can be adjusted through the removable elastic Velcro® Stabilo Strap
- The lengthwise elastic stretch zone in the elbow area ensures an ideal fit at all bending angles
- Moisture wicking fabric and antimicrobial finish wicks perspiration away, keeping the arm cool during activity and eliminating odour
- The top and bottom edges are made of extra soft snag proof material for extra comfort
- Compression class II (23-32 mmHg).

**SIZING:** Circumference 5cm below elbow joint

Size	Forearm Circumference
XS	8.7" - 9.4"
S	9.4" - 10.2"
M	10.2" - 11"
L	11" - 11.8"
XL	11.8" - 12.6"

**Indications:**

- Overstrain prophylaxis
- Sprains and strains
- Irritations caused by overstraining



## 816 Sportec Arm Compression Sleeve

**Features:**

- Special knitting technique that targets the forearm and bicep muscles to enhance strength
  - Reinforced around elbow joint to provide mild support for injury prevention, compressing tendons and ligaments
  - Medical grade compression of 20mmHg
  - Maximum pressure at the distal extremity and progressively decreasing to the proximal edge of the sleeve to improve circulation and remove lactic acid from muscles
- Sold in singles.

**SIZING:**  
(in cm)



Size	XS	S	M	L	XL	XXL
cX	19-26	22-29	25-32	28-36	33-40	37-44
cY	17-24	20-27	23-30	26-33	29-36	32-39
cZ	11-13	13-15	15-17	17-19	19-21	21-24



## 6011 Foam Padded Arm Sling/Shoulder Immobilizer

**Features:**

- Suspends arm and hand following injuries
- Soft padding over the shoulder and across the waist for comfortable wear
- Waist strap to stop upward motion of arm and shoulder rotation
- Can be used without waist strap for simple suspension

**SIZING:** Universal

### 5860 Ultra Cubital Tunnel Elbow Splint

**Features:**

- Immobilizes and positions elbow to reduce pain from cubital tunnel syndrome.
- Minimizes stretching of the Ulnar nerve.
- Cost-effective alternative to pre or post cast application.
- Made of a rigid black high density polyethylene shell with plush black terry cloth lining. Replacement liners also available.

Product #	Size	Description
5860SM	S/M	Fits most females and smaller framed males
5860ML	M/L	Fits most males and larger framed females
5860LXL	L/XL	Fits large males and XL framed females



### 531 SA Sling (Shoulder/Arm Immobilizer)

**Features:**

- Envelope is constructed of high quality black mesh for air circulation
- Padded foam shoulder strap is adjustable to suit any body type
- Shoulder immobilization waist strap is removable when only arm suspension is required
- Thumb loop provides additional support to the wrist

**Indications:**

- Cast support
- Shoulder separation
- Shoulder dislocation
- Shoulder/arm injuries

**SIZING:**

Size	Elbow to Fingertips	Waist
S	< 15"	< 32"
M	15" - 18"	< 40"
L	18" - 21"	< 48"

Alternative item for #5600 and #6014.



### 23191 Shoulder Abduction Pillow with Sling

**Indications:**

- Bankart repairs
- Rotator cuff repairs
- Shoulder arthroscopy, separations, and dislocations

**Features:**

- Easy-to-apply design with V-lock strapping system to prevent rolling and shifting
- Constructed with firm density foam covered in tricot material with tricot sling and exercise ball to stimulate circulation

**SIZING:**

Size	Waist Circ.	Elbow to Wrist
S	>33"	>8"
M	>38"	8" - 10"
L	>43"	10" - 12"
XL	>48"	10" - 14"



### 540 Shoulder Suspension Brace

**Indications:**

- Tendonitis and overuse injuries

**Features:**

- Adjustable strap to allow more influence over shoulder positioning
- Adjustable in the back and front to allow for a universal left and right fit
- Non-slip lining and adjustable opening at the bicep ensures the brace does not move from the desired position

**SIZING:**

Size	Bicep Circumference	Chest Circumference
S-M	10.5" - 12.5"	25" - 41"
L-XL	12.5" - 15.5"	41" - 54"





Product No.	Product	Page
-------------	---------	------

**Sacro-Iliac Supports**

902	Active SI Belt	D2
901	Universal Trochanter Belt	D2

**Pregnancy Supports**

904	Maternity SI Belt	D2
3090	Comfy Cradle Maternity Support	D2

**Lumbosacral Supports**

3250	Lumbosacral Back Support	D3
5549	Elastic Lumbo-Sacral Back Support	D3
6040	Lace Back Support	D3
580101-06	Aspen Lumbar	D3
6553	Pavis Ortho-E Back Brace	D4
5540	Dynamic LSO	D4
5543	Graphite Lite LSO	D4
5561	Edge LSO	D4

Product No.	Product	Page
-------------	---------	------

**Abdominal Binders**

1930 -1948	Surgical Binder	D5
6007	Abdominal Support (20cm)	D5
6009	Abdominal Support (32cm)	D5
104070	Bort Umbilical Hernia Support	D5
104090	Bort Abdominal Support w/ Ostomy	D5
6611	Pavis Rib Cage Binder	D6
6678	Pavis Wellness Ostomy Support Belt	D6
6685	Pavis Abdominal/Umbilical Hernia Support	D6

**Clavicle/Posture Supports**

6165	Front Closure Posture Corrector	D7
6001	Posture Back Support	D7
6164	EZ Adjust Clavicle Strap	D7

**TLSOs**

5545	Euro TLSO	D7
------	-----------	----



## 902 SI Strap

Designed to provide compression and stability to the sacro-iliac joints in a comfortable, low profile design.

**Features:**

- 3" wide strap features removable sacral pad for localized support
- Duo adjustable elastic straps allow for comfortable adjustment of focused support according to patient needs
- Strategically placed anti-migration pads to maintain brace position

**SIZING:**

Size	Hip Circumference
S	30"-36"
M	32"-38"
L	36"-44"
XL	40"-48"
XXL	46"-54"



## 901 Universal Trochanter Belt

Designed to gently compress the sacroiliac joint and hips. 3" wide elastic with Dura-Foam side panels and a hook and loop closure. Features adjustable elastic side pulls over the sacral area for added tension and stability.

**SIZING:** Universal (28" - 50" hips)



## 904 Padded SI Belt

**Indications:**

- SI joint pain relief
- Low back and hip pain relief during pregnancy and postpartum

**Features:**

- Comfortable compression and relief without bulk
- Lightweight, breathable material
- Conceals under clothing to firmly support pelvic structures
- No buckle design with detachable pad for extra support

**SIZING:**

Size	Hip Circumference
S/M	Up to 36"
L/XL	Up to 45"



## 3090 Comfy Cradle Maternity Lumbar Support

**This item has been discontinued, limited stock available. Please contact to confirm.**

The 8" pushback elastic pregnancy belt helps support the lower abdomen and back, offering relief from discomfort. The duo-adjustable, elastic side tension straps provide additional support, as does the optional mouldable insert. Use only under supervision and prescription of health care provider.

**SIZING:** Current measurement while pregnant

Size	Hip Circ.	Clothing Size
S/M	30" - 43"	4-12
L/XL	37" - 48"	12-20

### 3250 Lumbosacral Back Support

**Indications:**

- Chronic low back pain
- Low back strains/sprains
- Arthritic condition

**Features:**

- Generates support to the abdomen and lumbosacral areas
- 9" high elastic back panel has duo-tension straps that extend from the back and stretch to desired tension, adhering by Velcro
- Six flexible metal stays are enclosed to give added stability and prevent rolling

**SIZING:**

Size	Hip Circumference
S	22" - 28"
M	27" - 33"
L	32" - 38"
XL	37" - 43"
XXL	42" - 48"

Please note new sizing as of 2026.



### 5549 Elastic Lumbo-Sacral Back Support

**Features:**

- Ideal for stabilizing the lumbar spine, correcting posture, and providing pain relief
- Advanced German knitting technique with 4-way stretch, even compression, and breathability
- Anatomical knitting ensures a proper and comfortable fit to contour to the body
- Removable, injection-molded pad provides a massage effect, improves blood circulation to reduce muscle tension and pain

**SIZING:**

Size	Waist Circumference
S	81 - 96cm
M	94 - 109cm
L	107 - 122cm
XL	120 - 135cm



### 6040 Lace Back Support

Designed to help relieve lumbar area back pain caused by strains and muscle fatigue.

**Features:**

- Dynamic pulley system provides smooth and controlled compressive support
- Permits full range of movement while providing firm yet comfortable support - perfect for daily activities
- Side stretch panels help minimize bunching and provide easy for donning
- Latex and Neoprene free

**SIZING:**

Size	Waist Circumference
S/M	24" - 33"
L/XL	33" - 43"
XXL	43" - 50"



### Aspen Lumbar



Designed for active patients who require effective pain relief.

**Features:**

- Elegant low-profile brace for easy use under clothing, while still retaining a tailored, comfortable fit
- Aspen's patented 4 to 1 compression system allows independent adjustment of the top and bottom of the brace for customized fit
- Back height is 9", front height is 5.5"

**SIZING:**

Size	Waist Circumference	Product #
XS	53-69 cm	<b>580101</b>
S	66-79 cm	<b>580102</b>
M	79-91 cm	<b>580103</b>
L	91-107 cm	<b>580104</b>
XL	104-122 cm	<b>580105</b>
XXL	109-132 cm	<b>580106</b>





### 6553-28 Pavis Ortho-E Back Brace

**Features:**

- Semi-rigid textile in the paravertebral area is used to provide high level of column stability
- High quality elastic front straps effectively improve support
- Back stays are removable to adjust to patient's necessities and different rehab phases
- Innovative bacteriostatic fabric promotes hygiene, comfort, breathability, and prevents bad odours
- 28cm height, other height are available for special order, please contact us

**SIZING:**

Size	Pelvis Circumference
S	80 - 90 cm
M	90 - 100 cm
L	100 - 110 cm
XL	110 - 120 cm
XXL	120 - 130 cm
XXXL	> 130 cm



### 5540B Dynamic LSO

**Features:**

- Dual independent mechanical advantage pulley systems for precise upper and lower compression and maximum support
- Flexible pulley panels contour to the natural spinal curve
- Easy to use drawstring tabs for low hand strength clients
- Built-in abdominal panel for intra-abdominal pressure
- Removable pulley system makes washing easy
- Back height is 8", front height is 6.5"

**SIZING:**

Size	Waist Circumference
S	28-34"
M	34-40"
L	40-46"
XL	46-52"



### 5543 Graphite Lite LSO

**Indications:**

- Unstable lower spine, spasming muscles, and lower back pain

**Features:**

- Low-profile back panel and mesh fabric for cool comfort
- Graphite back panel contours to the shape of the spine to provide direct compression to the lower back
- EZ pull system has unique pull tabs for easy gripping, making it effortless to achieve effective compression
- Panel and lace system are removable for eash washing
- Back height is 10", front height is 5.5"

**SIZING:**

Size	Waist Circumference
S/M	24" - 33"
L/XL	40" - 50"



### 5561 Edge LSO

**Features:**

- Rigid spinal system to provide maximum stability
- Intracavitary pressure to unload intervertebral discs
- Back panel opening for wound inspection
- Comfortable moisture wicking padding
- Posterior height: 14.5", LSO Anterior height: 9", LSO LP anterior height: 7.5"

Product #	Description
<b>5561</b>	LSO
<b>5560</b>	LSO Low Profile (LP)

Size	Waist Circumference	Hip Circumference
S	30-36"	35-41"
M	34-40"	39-45"
L	38-44"	43-49"
XL	42-48"	47-53"



## Surgical Binder

**This item has been discontinued, limited stock available. Please contact to confirm.**

**Indications:**

- Pre or post abdominal surgery
- Mild abdominal hernias or muscle strains

**Features:**

- Provides compression and support
- 3 panel binder height is 9.5"
- 4 panel binder height is 12.25"
- For short-term use only

**SIZING:** Waist circumference

Product #	Description
<b>1930</b>	3 Panel Surgical Binder 22.5" - 45"
<b>1936</b>	3 Panel Surgical Binder 46" - 62"
<b>1940</b>	4 Panel Surgical Binder 30.5" - 45"
<b>1946</b>	4 Panel Surgical Binder 46" - 62"
<b>1948</b>	4 Panel Surgical Binder 62" - 78"



## 6007 Abdominal Binder (23cm)

**Indications:**

- Post abdominal surgery
- Pendulous abdomen

**Features:**

- Provides compression and support
- Durable design for long-term use
- 23cm high

**SIZING:**

Size	Waist Circ.
S	22" - 28"
M	27" - 33"
L	32" - 38"
XL	37" - 43"
XXL	42" - 48"



## 6009 Abdominal Binder (32cm)

The abdominal binder is designed to provide compression and support to the abdominal region after surgery. Its durable design also makes it ideal for long-term use with pendulous abdomens associated with obesity.

32cm high.

**SIZING:**

Size	Waist Circ.
S	22" - 38"
M	27" - 33"
L	32" - 38"
XL	37" - 43"
XXL	42" - 48"
XXXL	47" - 53"



## 104070 Bort Umbilical Hernia Support


**Features:**

- Silicone pad with plastic reinforcement
- Adhesive layer on the pad allows individual positioning
- Elastic fabric and two supporting splints
- Height of 13cm

Replacement silicone pads are available (**104080**)

**SIZING:**

Size	Waist Circumference
1	60-95 cm
2	95-125 cm
3	125-150 cm



## 104090 Bort Abdominal Support with Ostomy


**Features:**

- Reinforcing plastic panel designed to manage parastomal hernias
- Circular opening is adjustable in size and velour-lined panel is removable
- Height of 22cm
- Bag is not included

**SIZING:**

Size	Waist Circumference
S	70-90 cm
M	90-115 cm
L	115-150 cm





## 6611 Pavis Rib Cage Binder

**Indications:**

- Post chest trauma resulting in cracked or fractured ribs
- Post thoracic surgeries including heart surgeries

**Features:**

- Double adjustable crisscrossing frontal straps allow a correct, progressive, and comfortable closure
- Anatomically shaped contours around the armpit maximize comfort and compression to the upper chest
- Two adjustable straps stabilize the belt and prevent bending
- Bacteriostatic, hypoallergenic, and highly elastic fabric

**SIZING:**

Size	Chest Circumference
S	75 - 95 cm
M	95 - 115 cm
L	115 - 140 cm



## 6678 Pavis Wellness Ostomy Support Belt

**Features:**

- 9" wide ostomy belt with adjustable elastic compression to prevent hernias
- Colostomy hole has a 7cm diameter, can be adjusted up to 11 cm by cutting in the circular secured front panel
- 100% cotton fabric, combined with bacteriostatic and hypoallergenic fabric for maximum comfort
- Removable bag cover in elastic cotton attached to front

**SIZING:**

Size	Waist Circumference
S	75-90 cm
M	90-115 cm
L	115-135 cm



## 6685 Pavis Abdominal/Umbilical Hernia Support

Designed to contain the abdomen and the abdominal/umbilical hernia.

**Features:**

- Included pad can be easily repositioned as needed
- Can be used with pad before surgery or without after
- Wide velcro closure combined with woven elasticity allows tension adjustment to the desired level to reduce pain, improve efficiency, and providing stability
- Made with breathable cotton and elastic for highest comfort

**SIZING:**

Size	Hip Circumference
S	70 - 82 cm
M	82 - 94 cm
L	94 - 105 cm
XL	105 - 117 cm
XXL	117 - 130 cm

### 6165 Front Closure Posture Corrector

**Indications:**

- Poor posture
- Weak upper back muscles
- Pinched nerves
- Hernias of intervertebral disks
- Osteochondrosis
- Spondylosis of the spine


**SIZING:**

Size	Waist Circumference
XXS (Ped)	52-63 cm
XS (Youth)	63-69 cm
S	68-76 cm
M	76-85 cm
L	82-95 cm
XL	88-105 cm


**Features:**

- Retracts shoulders, provides support to clavicle
- Vest style is adjustable in the front, easy to don and doff
- Low profile design easily concealed under clothing

### 6001 Posture Back Support

**Indications:**

- Shoulder weakness
- Poor posture
- Clavicle injury rehabilitation


**SIZING:**

Size	Rib Circumference
S	20" - 28"
M	26" - 34"
L	32" - 40"
XL	38" - 46"


**Features:**

- Elastic shoulder straps connect to a soft cotton flannel back
- Front closure allows easy application
- Lightweight and discreet, does not have stays

### 6164 EZ Adjust Clavicle Strap

**Indications:**

- Pinched nerves
- Hernias of intervertebral disks
- Osteochondrosis
- Spondylosis of the spine

**Features:**

- Adjustable in the front
- Adjustable straps provide a comfortable fit for patients of all sizes
- Can be used for rehabilitation after injuries and spinal operations


**SIZING:** Universal


### 5545 Euro TLSO


The Euro TLSO provides rigid thoracic support to the lumbar spine in order to correct poor thoracic alignment. The rigid plastic spine stay is pre-shaped to fit the back and can be removed for custom shaping if desired. The advanced welding process creates a streamlined profile and a great fit. Shoulder pads provide extra comfort to the already soft fabric. Has a 5" front and 19" back.

**Indications:**

- Compression fractures
- Kyphosis
- Osteoporosis

**SIZING:**

Size	Waist Circumference
S	56-84 cm
M	81-109 cm
L	107-122 cm





Product No.	Product	Page
-------------	---------	------

### Hip/Groin Supports

201100	Bort StabiloHip Protector	E2
--------	---------------------------	----

### Hernia Supports

104800	Bort Hernia Belt	E2
6228	Hernia Support w/Pad	E2
6934	Pavis Prostatic Cushion	E2
6939	Pavis Coccyx Cushion	E2
6655	Wellness Boxer Hernia Support	E3
6656	Wellness Brief Hernia Support	E3
6648	Pavis Hernia Block Extra Strong with Pads	E3
6606	Pavis Replacement Hernia Support Pad	E3

### Ostomy Supports

6691	Ostomy Underwear for Men	E4
6692	Ostomy Underwear for Women	E4



## 201100 Bort StabiloHip Protector



Retail

This high quality hip protector is designed to protect the geriatric hip from fracture due to falling. Made from cotton and Lycra with a highly elastic waistband and leg borders. Bilateral protectors are made of shock absorbing visco-elastic soft foam and are easily removed.

Product #	Description
201100	Stabilo Hip Protector
201120	Replacement Pads
201200	Briefs Only
201250	Stabilo Hip Incontinence

### SIZING:

Size	Waist Circumference	Hip Circumference
S	27" - 29"	35"
M	27" - 31"	37"
L	31" - 32"	38"
XL	32" - 34"	40"
XXL	34"+	42"+



## 104800 Bort Hernia Belt



This high quality Bort Hernia Belt has a flexible anterior panel to nicely conform to the abdominal area. All of the straps are made from high quality, strong elastic to provide excellent compression. The leg straps hold the hernia pads in place for optimal positioning of the hernia support pads. The pads are removable for single sided use and attach with Velcro® to be placed in any position. Indicated for reducible inguinal hernia.

### SIZING:

Size	Hip Circ.
S	up to 89 cm
M	89-104 cm
L	104-119 cm



## 6228 Hernia Support with Pads

Provides firm support and comfort, and fits under clothing. The elastic strap is 2 inches wide. The flexible front panel fits body contours and the leg straps are made of soft, knitted elastic. Has right and left moulded pads for support. Pads are removable for single sided use and attach with Velcro®.

### SIZING:

Size	Hip Circ.
S	30" - 35"
M	35" - 41"
L	41" - 46"
XL	46" - 51"



## 6934 Pavis Prostatic Cushion

Made of viscoelastic foam with pressure relief points strategically positioned to adapt perfectly to the user's anatomy. A hollow central cavity promotes ventilation from the perianal to the genital area, avoiding any contact and preventing pressure. Indicated for people with prostatic or hemorrhoidal conditions, cases of pelvic or anorectal surgeries with low sperm counts, pregnant women, mothers in postpartum recovery, or people who remain seated for long periods of time. Has a waterproof cover.

**Dimensions (cm):** 42 x 36 x 8.5



## 6939 Pavis Coccyx Cushion

The "U" shaped curve allows the decompression and ventilation of the anorectal and perianal areas. Indicated for people with hemorrhoidal conditions, subject to pelvic and anorectal surgeries, or for those who remain seated for long periods of time. Has a velvet cover.

**Dimensions (cm):** 42 x 36 x 8.5

### 6655 Wellness Boxer Hernia Support

The Boxer Hernia Support integrates the design of a traditional hernia support with high-strength elastic and removable hernia support pads into a stylish and comfortable pair of boxer shorts. The hernia pads are anatomically shaped to reduce the hernia for pain relief. High-tech athletic materials are used to provide breathability, antibacterial action and high durability for full-time wear and washing. Nanotechnology yarn ensures the bacteriostatic quality remains unchanged for more than fifty washing cycles. The boxer style is very effective for active use during sport activities.

**SIZING:**

Size	Pelvis Circumference
S	71-82 cm
M	82-94 cm
L	94-105 cm
XL	105-117 cm
XXL	117-130 cm



### 6656 Wellness Brief Hernia Support

This brief-style hernia support integrates the design of a traditional hernia support with high strength elastic and removable hernia support pads into a stylish and comfortable pair of briefs. The hernia pads are anatomically shaped to reduce the hernia for pain relief. High-tech athletic materials are used to provide breathability, antibacterial action and high durability for full time wear and washing. Nanotechnology yarn ensures the bacteriostatic quality remains unchanged for more than fifty washing cycles. The boxer style is very effective for active use during sport activities.

**SIZING:**

Size	Pelvis Circumference
S	71-82 cm
M	82-94 cm
L	94-105 cm
XL	105-117 cm
XXL	117-130 cm



### 6648 Pavis Hernia Block Extra Strong with Pads

Similar to the popular Pavis Hernia Briefs, we offer these extra strong hernia supports for the hernia that is more difficult to hold. A set of left and right hernia pads are included to contain reducible inguinal hernias. Briefs are also recommended after hernia surgery to prevent abdominal relaxation, and in that case should be worn without the pads. This model is high on the abdomen and lower on the thigh to provide optimum support and containment. Maximum functionality and wearability is guaranteed by the use of different types of fabric within the product, all characterized by different elastic forces. Unique breathable 100% cotton elastic fabric in contact with the skin, made exclusively by Pavis, provides excellent comfort and is non-allergenic.

**SIZING:**

Size	Pelvis Circumference
65	71 - 77 cm
70	77 - 82 cm
75	82 - 88 cm
80	88 - 94 cm
85	94 - 100 cm
90	100 - 105 cm
95	105 - 110 cm
100	110 - 117 cm
105	117 - 123 cm
110	123 - 130 cm



### 6606 Pavis Replacement Hernia Support Pad

Replacement pads for the popular Pavis hernia supports are anatomically shaped to contain a hernia in combination with the Pavis hernia supports. They are molded from firm foam with a velour cover. Specify left or right.





## 6691 AirGeniX Ostomy Underwear for Men

**Features:**

- The first Pavis boxers designed for stoma containment
- High-waisted boxer-type underwear, with a modern design
- Two internal pockets, which can be used to contain the ostomy bag, keeping it hidden and supported inside the boxers
- Front opening for maximum effectiveness and comfort
- Made with the innovative AirGeniX antibacterial fabric, ensuring maximum hygiene, comfort, freshness, breathability, and preventing odours
- To be worn directly in contact with the skin as an undergarment (ostomy bag not included)

**SIZING:**

Size	Waist Circumference
S	76-84 cm
M	84-92 cm
L	92-100 cm
XL	100-108 cm
XXL	108-116 cm



## 6692 AirGeniX Ostomy Underwear for Women

**Features:**

- The first Pavis briefs designed for stoma containment
- High-waisted underwear, with a modern design
- Two internal pockets, which can be used to contain the ostomy bag, keeping it hidden and supported inside the briefs
- Made with the innovative AirGeniX antibacterial fabric, ensuring maximum hygiene, comfort, freshness, breathability, and preventing odours
- To be worn directly in contact with the skin as an undergarment (ostomy bag not included)

**SIZING:**

Size	Waist Circumference
XXS	58-66 cm
XS	66-74 cm
S	74-83 cm
M	83-92 cm
L	92-102 cm
XL	102-112 cm
XXL	112-122 cm
XXXL	122-134 cm
XXXXL	134-146 cm

Product No.	Product	Page
-------------	---------	------

### Wrist Compression

1361	Slip-On Wrist Compression	F2
5591	3D Elastic Wrist Support	F2

### Wrist Supports with Stays

87	Arthritic Anti-Ulnar Deviation Splint	F2
5593	Dynamic Wrist Stabilizer	F2
93N	Neoprene Wrist Lacer	F3
3193	CT Lacer	F3
3793	Airflex CT Lacer	F3
97A	Contoured Wrist Stabilizer	F4
95A	Wrist Stabilizer	F4
3795	Airflex Wrist Stabilizer	F4
6133	Pavis P-Fix Wrist Splint	F4
168	Exoform Carpal Tunnel Brace	F5
133300	Bort Wrist Extension Support	F5
112770	Bort DigiSoft Finger Brace (Boxer Splint)	F5

### Wrist and Thumb Supports

1362	Elastic Wrist Thumb Support	F5
69A	Neoprene Wrist Support with Thumb Extension	F6
69B	Thumb Stabilizer	F6
67	Wrist/Thumb Stabilizer	F6
3767	Airflex Thumb Support	F6
6635	Pavis Thumb Brace	F7
5567	Dynamic Thumb Stabilizer	F7
112080	Bort SellaFix MCP Stabilizer	F7
112040	Bort SellaFix CMCP Stabilizer	F7
105400	Bort SellaDur Stabilizer	F8
93NT	Neoprene Wrist Thumb Lacer	F8
3193ST	Coolcel Wrist Short Thumb Lacer	F8
3193T	Wrist Thumb Lacer	F8
6637	Pavis Wrist Thumb Lacer	F8

Product No.	Product	Page
-------------	---------	------

### Arthritic Gloves and Supports

ST-6838	Swede-O Thermal Arthritis Gloves	F9
5189	Thermoskin Full Finger Arthritis Gloves	F9
ST-6839	Swede-O Thermal Carpal Tunnel Glove	F9

### Finger Supports

9110	Resting Hand Orthosis	F10
564-568	Aluminum Padded Finger Splints	F10
572-578	Aluminum Finger Splints	F10

### Silicone Products

2600	DynaGel™ Digital Tube	F11
2626	DynaGel™ Digit Tube	F11
2601	DynaGel™ Digital Cap	F11
2633	DynaGel™ Digital Cap	F11



## 1361 Slip On Wrist Compression

**This item has been discontinued, limited stock available. Please contact to confirm.**

Beige knitted cotton elastic wristlet slips on for easy application. 6" long with thumb hole. Fits left or right hand.

SIZING:	Size	Wrist Circumference
	S	5.5" - 6.5"
	M	6.75" - 7.5"
	L	7.75" - 8.5"
	XL	8.75" - 10"



## 5591 3D Elastic Wrist Support

- Manufactured with 3D knitting technology to provide superior fit and comfort
- Provides medical compression of 20-30mmHg
- Features an elastic wrist stabilizing strap that provides added support to the active wrist for pain and tendon strains
- Knitted construction is breathable for cool comfort

SIZING:	Size	Wrist Circ.
	S	5" - 6"
	M	6" - 7"
	L	7" - 8"
	XL	8" - 9"
	XXL	9" - 10"



## 87 Arthritic Anti-Ulnar Deviation Splint

- Made with Neoprene to keep the arm and hand warm and to stimulate circulation
- Removable aluminum stays come pre-formed to support the wrist at the palm and on the medial ulnar aspect
- For the best support and ease of application, straps include a pull-pack with D-ring at the wrist and pull-over Velcro at the proximal arm opening
- Finger positioning straps pull dorso-radially and attach with Velcro on the palm
- Glove style wrist wrap provides easy donning, and palm extension tab provides easy doffing

**SIZING:** See chart below, specify left or right



## 5593 Dynamic Wrist Stabilizer

**Features:**

- Lace closure system for optimal fit and easy application
- Breathable foam material for cool comfort
- Palm and dorsal stays for good support
- Reversible stays to fit either left or right
- Only 3 sizes to fit a wide range of wrist sizes

**SIZING:**

Size	Wrist Circumference
S	13-17 cm
M	17-21 cm
L	21-23 cm

SIZING:	Size	Wrist Circumference	Size	Wrist Circumference
	XS	4" - 5"	L	7" - 8"
	S	5" - 6"	XL	8" - 9"
	M	6" - 7"	XXL	9" - 10"

### 93N Neoprene Wrist Lacer

The Neoprene Carpal Tunnel Wrist Lacer is the Neoprene version of our popular Active CT Lacer. It is designed to stabilize weak, unstable wrists and generate warmth to stimulate circulation. The easy-to-use Velcro® closure lacing system allows for easy donning and doffing – the brace can be adjusted with one simple motion.

**SIZING:** See chart on opposite page, specify left or right.

Product #	Description
93N	6.5" Length
93LN	8" Length



### 3193 Coolcel CT Lacer

The CT Lacer is the ideal choice for most patients. With an easy-to-use Velcro® closure lacing system, the brace may be easily donned with one simple motion. The lacing system is designed so that the proximal and distal straps may be adjusted independently of one another, allowing for a proper fit for patients with unique hand and wrist sizes.

This brace is designed to stabilize weak, unstable wrists. The unique Ortho Active Contour dorsal and palmar stays exhibited below are incorporated into this brace. These aluminum stays hold the wrist in a neutral, cock-up position for proper healing. They are removable and adjustable, and provide a low profile, intimate, and comfortable fit. The brace material is cut proximal to the palm crease to allow full hand motion while maintaining immobilization of the wrist. Comes in 6.5" length, also available in 8" length (3194). Specify left or right.

#### Ortho Active Contour Stays

Dorsal stay with contour relief for styloid process



Palmar stay with thumb contour



#### Indications:

- Carpal tunnel syndrome, tendonitis
- Post-cast
- Moderate strains or sprains

**SIZING:** See chart on opposite page, specify left or right

Product #	Description
3193	6.5" Coolcel Lacer
3194	8" Coolcel Lacer



### 3793 Airflex Carpal Tunnel Wrist Lacer

The Active Airflex CT Lacer is the ultimate in comfort, fit, and ease of use. It is designed to provide compression and support to the affected wrist. The stays are anatomically contoured to provide a low-profile fit that does not impede use of the hand or fingers. Lacer design is easy to apply.

**SIZING:** See chart on opposite page specify left or right.



## 97A Neoprene Contoured Wrist Stabilizer

Provides stability for weak/unstable wrists. Contour distal and palmar aluminum stays, that are removable and adjustable, hold the wrist in a neutral position for healing. They provide a low profile and comfortable fit. The brace material is cut proximal to the palm crease to allow full hand motion while maintaining immobilization of the wrist. Application is easy with the aid of a d-ring. Available in Neoprene or Coolcel (ideal for sport use or patients who are prone to irritation). 6.5" length.

**SIZING:** See chart below, specify left or right

**Indications:**

- Carpal tunnel syndrome, tendonitis
- Post-cast
- Arthritic wrist
- Moderate strains or sprains

Product #	Description
97A	6.5" Neoprene
3197A	6.5" Coolcel



## 95A Wrist Stabilizer

Improved, industrial strength, 8" wrist stabilizer with a D-ring tension strap to stabilize weak, unstable wrists and generate warmth (Neoprene) to stimulate circulation. Two mouldable aluminum stays hold the wrist in a neutral, cock-up position for proper healing. Thumb hole is adjustable.

**Indications:**

- Carpal tunnel syndrome, tendonitis
- Post-cast
- Arthritic wrist
- Moderate strains or sprains

**SIZING:** See chart below, specify left or right.

Product #	Description
95A	8" Neoprene
3195A	8" Coolcel



## 3795 Airflex Industrial Wrist Brace

The Ortho Active Airflex Double Stay Industrial Wrist Brace is designed to provide stability to the injured wrist. The stays may be removed and molded to provide support as needed. The 8" length and wide straps are ideal for the working hand and wrist.

**SIZING:** See chart below, specify left or right.



## 6133 Pavis P-Fix Wrist Splint

**Indications:**

- Trauma/sprains, injury of flexor and extensor tendons
- Tendinitis, osteoarthritis, arthritis, carpal tunnel

**Features:**

- 3 anatomically shaped, mouldable and removable stays
- 3 straps to adjust closure according to patient measurements
- Totally open brace for easy wearing with pain or after operation
- Antibacterial AirGeniX textile to improve hygiene, comfort, and breathability
- 22cm length, for left or right

**SIZING:**

Size	Wrist Circ.
Small	15 - 17cm
Medium	17 - 21cm
Large	21 - 25cm

SIZING:	Size	Wrist Circumference	Size	Wrist Circumference
		XS	4" - 5"	L
	S	5" - 6"	XL	8" - 9"
	M	6" - 7"	XXL	9" - 10"

### 168 Exoform Carpal Tunnel Brace

**This item has been discontinued, limited stock available. Please contact to confirm.**

The lightweight, low profile design promotes good patient compliance, the contoured exoskeleton ensures an anatomically correct fit, and provides smooth controlled resistance to extreme range of motion. It is pre-shaped into a neutral position for comfortable, limited motion/flexibility. Provides radial and ulnar styloid relief, and more room at the base of the thumb, palmar crease and 5th MP.

**Indications:**

- Carpal tunnel syndrome
- Tendonitis
- Moderate strains or sprains

**SIZING:** Specify left or right

Size	Wrist Circumference	Palm Circumference
S	5.25" - 6.5"	6.75" - 8"
M	6.5" - 7.75"	8" - 9.25"
L	7.75" - 10"	9.25" - 10.5"



### 133300 Bort Wrist Extension Support



This elastic support is for immobilization in the case of hand and finger injuries, and for tenosynovitis of the fingers.

**Instructions for use:** Bend aluminium splint to the desired form before putting it on. First close fastening between thumb and forefinger, then place Velcro® fastening over the fingers. Has three Velcro® straps.

**Indications:**

- Tenosynovitis
- Mild contractures

**SIZING:** Specify left or right

Size	Wrist Circ.
S	15 - 17 cm
M	17 - 19 cm
L	19 - 21 cm
XL	21 - 23 cm



### 112770 Bort DigiSoft® Finger Brace



The Bort DigiSoft® Finger Brace provides immobilization of the interphalangeal joints for either one or two fingers. It is made of adjustable aluminum, covered with soft padding for maximum comfort. The Velcro fasteners on the metacarpals and fingers can be adjusted as required.

**Indications:**

- Post-operative
- Post-traumatic
- Positioning

**SIZING:** Length from wrist to tip of longest injured finger

Size	Hand Length
1	13-16 cm
2	16-19 cm
3	19-22 cm



### 1362 Elastic Wrist Thumb Support

The Ortho Active Elastic Wrist Thumb Support is designed to provide mild compression to the wrist. It features durable, knitted cotton elastic, and it fits either left or right. This support is similar to our popular 1361 Elastic Wrist Support with the addition of a thumb support.

**SIZING:**

Size	Wrist Circ.
S	13-15 cm
M	15-17 cm
L	17-19 cm
XL	19-21 cm
XXL	21-23 cm





## 69A Universal Size Thumb Support

Same as 69 but with a soft thumb extension for added support. 3mm Neoprene.

**SIZING:** Universal

Product #	Description
69A	Neoprene w/ thumb extension



## 69B Thumb Stabilizer

Similar to the 69A with a metal thumb stay to limit motion of the thumb.

Product #	Description
69B	Thumb Stabilizer
69BU	Universal

**SIZING:**

Size	Wrist Circumference
S	5" - 6"
M	6" - 7"
L	7" - 8"



## 67 Wrist/Thumb Stabilizer

The thumb support was designed for easy application and moderate support of the carpal metacarpal and metacarpal phalangeal joints to reduce pain and promote healing. A pre-formed moldable stay is incorporated and the support is made of soft padded neoprene for added comfort. The user can simply pull it over the hand and tighten the wrist strap to support the joints of the thumb while still being able to use the distal thumb joint for everyday tasks.

**SIZING:** Specify left or right

Size	Wrist Circ.
S	5" - 6"
M	6" - 7"
L	7" - 8"
XL	8" - 9"



## 3767 Airflex Thumb Support

The Airflex Thumb Lacer immobilizes the thumb with a pre-formed, moldable, aluminum stay. This comfortable support stabilizes the CMC and MCP joints. It is easy to apply with Velcro/D-ring closure and is adjustable for each individual with a Velcro closure over the thumb.

**SIZING:** Specify left or right

Size	Wrist Circ.
XS	4" - 5"
S	5" - 6"
M	6" - 7"
L	7" - 8"
XL	8" - 9"

### 6635 Pavis Thumb Brace

The Pavis Thumb Immobilizer incorporates the Pavis Rizofix patented system with adjustable Velcro closure, allowing the thumb to be immobilized in complete comfort. The dorsal thumb area is open for scar inspections or bandage care. It can be opened completely, making it easy and painless to apply and infinitely adjustable, ensuring a comfortable position. 100% soft cotton velour in contact with the skin, making this brace comfortable, breathable and hygienic. Fits left or right.



### 5567 Dynamic Thumb Stabilizer

Features:

- Lace closure system for optimal fit and easy application
- Breathable foam material for cool comfort
- Removable mouldable thumb stay for good support
- Fits either left or right hand

**SIZING:** Specify left or right

Size	Wrist Circumference
S	13-17 cm
M	17-21 cm
L	21-23 cm



### 112080 Bort SellaFix MCP Stabilizer



A rigid plastic orthosis for stabilizing the MCP joint while allowing full freedom of movement at the wrist in the case of sprains and sprains. Lightweight and simple to mould and cut.

**Indications:**

- Sprains and strains
- Skier's thumb
- Gamekeeper's thumb

**SIZING:** Specify left or right

Size	Wrist Circumference
XS	up to 15 cm
S	15 - 17 cm
M	17 - 19 cm
L	19 - 21 cm
XL	21 - 23 cm



### 112040 Bort SellaFix CMCP Stabilizer



A rigid plastic orthosis for stabilizing the CMC joint and the MCP joint. No impairment of grip function and fine motor activity. Simple to mould and cut, lightweight.

**Indications:**

- Ulnar ligament lesions (skier's thumb)
- Post scaphoid fracture
- Post surgery
- Degenerative joint disease in thumb

**SIZING:** Specify left or right

Size	Wrist Circumference
XS	up to 15 cm
S	15 - 17 cm
M	17 - 19 cm
L	19 - 21 cm
XL	21 - 23 cm





### 105400 Bort SellaDur Stabilizer



An economical version of the SellaFix CMCP stabilizer. This thumb saddle joint orthosis with perforation is comfortable to wear. The special high-formed edge avoids pressure spots in the critical thumb-forefinger area. Immobilizes the CMC and MCP joints. See indications for SellaFix CMCP **112040**. Comes with a black undersleeve.

**SIZING:** Specify left or right

Size	Wrist Circumference
Small	15 - 17 cm
Medium	17 - 19 cm
Large	19 - 21 cm

### Wrist Thumb Lacer



This support from Ortho Active provides immobilization of the wrist and thumb, while allowing mobility of the phalangeal joint. Features a removable, mouldable palmar stay and a dorsal stabilizer at the wrist. The thumb is immobilized by a pre-shaped, adjustable aluminum splint. With an easy-to-use Velcro® closure lacing system, the brace may be donned with one simple motion. Comes in both Neoprene and Coolcel, Ortho Active's lightweight, latex-free, Neoprene alternative.

**Indications:**

- DeQuervain's Tenosynovitis
- Post cast after scaphoid fracture
- Thumb sprains, strains, or tendonitis
- Wrist conditions requiring thumb immobilization

**SIZING:** Specify left or right

Size	Wrist Circ.
XS	4" - 5"
S	5" - 6"
M	6" - 7"
L	7" - 8"
XL	8" - 9"
XXL	9" - 10"

Product #	Description
<b>93NT</b>	Neoprene Wrist Thumb Lacer
<b>3193ST*</b>	Coolcel Wrist Thumb Lacer



### 3193T Coolcel Wrist Thumb Lacer

This support from Ortho Active is similar to the **3193** CT Lacer but also provides immobilization of the thumb. Features a removable, mouldable palmar stay and a dorsal stabilizer at the wrist. The thumb is immobilized by a pre-shaped, adjustable aluminum splint. With an easy-to-use Velcro® closure lacing system, the brace may be donned with one simple motion. It is made with Coolcel, Ortho Active's lightweight, latex-free, Neoprene alternative.

**Indications:**

- DeQuervain's Tenosynovitis
- Post cast after scaphoid fracture
- Thumb sprains, strains, or tendonitis
- Wrist conditions requiring thumb immobilization

**SIZING:** See chart above for #3193ST



### 6637 Pavis Wrist Thumb Lacer

The three Velcro closures with adjustable anchor point for the laces make this an easy-to-fit wrist/thumb brace. Four mouldable, removable, pre-contoured stays provide excellent and comfortable support for the wrist and thumb. The lace/Velcro closure system is easy to adjust and tighten. Specify Left or Right.

**SIZING:** Specify left or right

Size	Wrist Circumference
Small	10-14 cm
Regular	14-20 cm
Extra	20-26 cm

### ST-6838 Swede-O Thermal Arthritis Gloves

- Allows skin to ventilate and remain well oxygenated for extended periods
- Increases blood flow, facilitating the removal of cellular debris
- Provides light but firm compression to counteract tissue swelling
- Promotes optimal muscle function
- Increases elasticity and reactivity of the muscle, so after its application the collagen fibres become more extensible and capable of rehabilitation exercises
- Provides temporary relief from pain/soreness associated with sports injuries, arthritis and RSI

**Comes in pairs**

**SIZING:** Circumference of the hand around the knuckle

Size	MCP Circ.
XS	6"-6.75"
S	7"-7.75"
M	8"-9"
L	9.25"-10.5"
XL	10.75 -11.5"
XXL	11.75"+



### 5189 Thermoskin Full Finger Arthritis Gloves



**This item has been discontinued, limited stock available. Please contact to confirm.**

Provides heat therapy and pain relief to the entire finger. Thermoskin is clinically proven to increase blood flow and skin temperature. Provide light compression and support to reduce soreness and swelling, providing natural relief.

Useful when doing household chores and gardening, and they are also helpful to anyone requiring gentle heat therapy to the joints of the hand. The palm area of the gloves are textured with Griptex to provide additional grip and the smooth outer layer is smooth nylon for appearance and ideal for when worn overnight to prevent catching on sheets while sleeping. There is a Velcro wrist closure for easy application and better fit. There is a right-hand and left hand glove in every pack.

**SIZING:** Circumference of the hand around the knuckle

Size	MCP Circ.
S	7"-7.75"
M	8"-9"
L	9.25"-10.5"
XXL	10.75" - 11"



### ST-6839 Swede-O Thermal Carpal Tunnel Glove

- Similar to the **ST-6838** but has a flexible palmar stay
- Helps prevent wrist movements that can cause Carpal Tunnel Syndrome
- Comfortable, firm support for weak, injured or aching wrists
- Provides relief of Carpal Tunnel Syndrome
- Allows for full use of fingers
- Added benefit of heat therapy
- Relieves pain and inflammation
- Anatomically conforming fit

**Comes in singles**

**SIZING:** Circumference of the hand around the knuckle, specify left or right

Size	MCP Circ.
S	7"-7.75"
M	8"-9"
L	9.25"-10.5"
XL	10.75"-11.5"




### 9110 Resting Hand Orthosis



This splint offers a functional resting hand position following stroke, trauma or injury, and can also be used in the treatment of arthritis, burns, and early wrist and finger contractures. Constructed of a lightweight, durable Kydex plastic, the Resting Hand Orthosis can be heat moulded at a low temperature to meet individual patient needs. The washable, breathable foam liner eliminates potential pressure points while wicking moisture away, keeping the skin dry and comfortable. Specify left or right.

#### SIZING:

Size	Palm Width	Hand Length
Small	up to 3"	up to 6"
Medium	3" - 3.5"	6.5" - 8"
Large	3.5" - 4"	8" - 9"




### Aluminum Finger Splint

This pre-cut, smooth, trimmed finger splint has foam padded aluminium which is bonded to unique, cushioned, self adhering straps which cradle the finger securely. Easy to use while providing the utmost comfort and stability.

Product #	Size	Length
<b>572</b>	S	4.5" Length
<b>574</b>	M	6.25" Length
<b>576</b>	L	7.5" Length
<b>578</b>	Universal	2.25" Length

### Aluminum Padded Finger Splints



Product #	Description	Measurement
<b>564</b>	Finger Foldover Splint	S ( 1.5"), M ( 2"), L ( 3.5")
<b>565</b>	Finger Protector Splint	S ( 2"), M ( 3"), L ( 3.5")
<b>566</b>	Finger Gutter Splint	S ( 1.5"), M ( 3"), L ( 4.5"), XL ( 6")
<b>567</b>	Finger Frog Splint	S ( 2.3"), M ( 3"), L ( 3.5")
<b>568</b>	Finger Baseball Splint	S-L

## 2600 DynaGel™ Digital Tube

- Mineral-oil-infused antibacterial gel provides padding, relieves pressure and eliminates friction. Medical grade mineral oil is slowly released to soften and moisturize the skin. Skin softening action will help to relieve scar tissue.
- Hyper-elastic fabric is conformable and comfortable for any shape.
- Long lasting gel will not wear out like foam products.

**2600-5** - Five 15cm tubes per bag

**R2600-4** - Four 2cm tubes per retail box

Size	Diameter
S	0.5 cm
M	1 cm
L	1.5 cm



## 2626 DynaGel™ Digit Tube

- Antibacterial gel provides padding, relieves pressure and eliminates friction.
- This soft, extra-durable gel does not need a fabric backing; it stretches easily and provides less compression than our regular ribbed tube. Long lasting gel will not wear out like foam products.
- Easily cut to fit required length.

**R2626-2** - Two 3cm tubes per retail box

Size	Diameter
S	0.5 cm
M	0.8 cm
L	1 cm



## 2601 DynaGel™ Digital Cap

- Mineral-oil-infused antibacterial gel provides padding, relieves pressure and eliminates friction. Medical grade mineral oil is slowly released to soften and moisturize the skin. Skin softening action will help to relieve scar tissue.
- Hyper-elastic fabric is conformable and comfortable for any shape.
- Long lasting gel will not wear out like foam products.
- Easily cut to fit required length.

**2601-10** - Ten 5cm caps per bag

**R2601-2** - Two 5cm caps per retail box

Size	Diameter
S	0.6 cm
M	0.8 cm
L	1 cm



## 2633 DynaGel™ Digital Cap

**This item has been discontinued, limited stock available. Please contact to confirm.**

- Antibacterial gel provides padding, relieves pressure and eliminates friction. This soft, extra-durable gel digital cap does not need a fabric backing, it stretches easily, and provides less compression than our regular ribbed cap. Long lasting gel will not wear out like foam products.
- Easily cut to fit required length.

**2633-10** - Ten 4.5cm caps per bag

**R2633-2** - Two 4.5cm caps per retail box

Size	Diameter
S	0.5 cm
M	1 cm
L	1.5 cm





Product No.	Product	Page
-------------	---------	------

### Ortho Active Knee Supports

32	Neoprene Hinged Knee Brace	G2
39	Neoprene Multi-Stop Hinged Knee Brace	G2
36	Spiral Knee Support	G3
34	Neoprene Jumpers Knee Brace	G3
29	Knee Brace Undersleeve	G3

### Patellar Strap

3126	Coolcel Infrapatellar Strap	G4
6026	Patella Tendon Strap	G4

### Knee Compression Sleeves

3511	Elastic Knee Support	G4
3512	Charcoal Elastic Knee Support	G4
5530	3D Elastic Knee Stabilizer	G4
3611	Slip-On Knee Compression	G5
820	Sportec Compression Sleeve	G5
3531	Sportec Patella Compression Knee Support	G5

Product No.	Product	Page
-------------	---------	------

### Knee Supports

5537	Dynamic Knee Support	G6
3623	Spiral Elastic Knee Support	G6
648101-648106	Aspen ROM Knee	G6
651021-651030	Aspen OA Knee Wrap	G6
150140	Bort Stabulo Knee Support with Hinge	G7
114420	Bort Select Stabilogen Knee Support	G7
114700	Bort Asymmetric Knee Support	G7
6640	Patella Stabilizer	G8
6642	Multi-Stability Lace Knee Brace	G8
6622	Hinge Knee Wrap	G8
6653	ACL-PCL/MCL-LCL Stabilizer Knee Brace	G8
5438	Post-Op Knee Brace	G9
122-125	Universal 3-Panel Knee Immobilizer	G9
11-4	Glyde OA	G9
06-5	TRUalign OA	G9
651011-651012	Aspen OA Knee+	G10
652001-652009	Aspen CI Knee	G10
651111	Suspension Wrap Strap	G10

### Sportec Sleeves

811	Sportec Calf Compression Sleeve	G11
805	Sportec Ankle Compression Sleeve	G11
815	Sportec Leg Compression Sleeve	G11



## 32 Hinged Knee Brace

The Ortho Active Hinged Knee Brace is designed to provide medial and lateral support to the knee as well as stabilize the patella and provide compression and warmth. The hinges may be removed and moulded to the leg to provide a custom fit. Available with an open-top or wrap-around style for easier brace application. Material options include Neoprene or Coolcel.

**Features:**

- Aluminum polycentric hinges with 10 degree hyperextension stops mimic anatomical knee motion
- Inferior/superior compression straps for extra support
- Removable patellar buttress assists in patellar tracking (not available in **32W** or **3132W**)

Product #	Description
<b>32P/3</b>	3mm Neoprene with Felt Buttress
<b>32PO/3</b>	3mm Neoprene with Felt Buttress - Open Top
<b>32D/3</b>	3mm Neoprene with Donut Buttress
<b>32W/3</b>	3mm Black Neoprene - Wrap Around (No Buttress)
<b>3132P</b>	Coolcel Hinged with Felt Buttress
<b>3132D</b>	Coolcel Hinged with Donut Buttress
<b>3132W</b>	Coolcel Hinged Wrap

**SIZING:** See chart below, also available in XXXL (18"+)



**Indications:**

- Medial/lateral collateral ligament instabilities
- Anterior/posterior cruciate ligament instabilities
- Minor meniscal injuries

## 39 Multi-Stop Hinged Knee Brace

The Ortho Active Multi-Stop Hinged Knee Brace is designed to provide medial and lateral support to the knee. It also stabilizes the patella, and provides compression and warmth. The hinges have adjustable flexion and extension stops that may be adjusted in 15° increments. Four Velcro® (2 inferior, 2 superior) straps allow for variable compression and work with the hinge to limit motion. Available in open-top and wrap-around styles for easy application. Material options include Neoprene or Coolcel.



**Indications:**

- Medial/lateral collateral ligament instabilities
- Anterior/posterior cruciate ligament instability
- Minor meniscal injuries

**SIZING:** See chart below. Also available in XXXL (18"+)

Product #	Description
<b>39P/3</b>	3mm Neoprene with Felt Buttress
<b>39PO/3</b>	3mm Neoprene with Felt Buttress - Open Top
<b>39W</b>	3mm Neoprene - Wrap Around
<b>3139P</b>	Coolcel Multi-Stop -With Felt Buttress
<b>3139W</b>	Coolcel Multi-Stop - Wrap Around

**Neoprene is available in black**

**SIZING:**

**Neoprene**

Always measure around the knee centre

XS	13" - 14"	L	16" - 17"
S	14" - 15"	XL	17" - 18"
M	15" - 16"	XXL	18" - 19"

**Coolcel**

Always measure around the knee centre

XS	12" - 13"	L	15" - 16"
S	13" - 14"	XL	16" - 17"
M	14" - 15"	XXL	17" - 18"

## 36 Spiral Knee Support

**Features:**

- Patella hole relieves pressure on the patella
- Three spring stays at medial/lateral aspects provide extra stability for weak ligamentous structures in the knee
- Provides warmth and stability, increasing circulation
- Allows function without causing further trauma

**SIZING:** See chart on opposite page

Product #	Description
36	3mm Neoprene Spiral Knee
3136	Coolcel



## 34 Jumpers Knee Brace

The Ortho Active Jumper's Knee Brace is designed to provide compression, warmth, and support to the knee. The inferior and superior straps reduce tension on the patellar tendons, thus reducing the pain associated with repetitive loading of the knee joint. Options include open-top and wrap-around styles for easy application. Available in Neoprene or Coolcel, Ortho Active's lightweight, latex-free Neoprene alternative. Specify left or right for **34J**.

**Features:**

- Inferior/superior compression straps for extra support
- Removable patellar buttress assist in patellar tracking and can be positioned for the needs of the user (not available in 34W or 3134W)
- Spiral stays prevent rolling and give medial/lateral support

**Indications:**

- Jumper's Knee (Patellar Tendonitis)
- Runner's Knee (Patellofemoral Pain Syndrome)
- Mild knee ligament sprains
- Arthritic knee

**SIZING:** Knee centre circumference  
See chart on opposite page

Product #	Description
34P/3	3mm Neoprene with Felt Buttress
34PO/3	3mm Neoprene with Felt Buttress-Open Top
34D/3	3mm Neoprene with Donut Buttress
34W/3	3mm Neoprene - Wrap Around (no buttress)
34J/3	3mm Neoprene with Lateral J
3134J	Coolcel Patella Lateral J Knee
3134P	Coolcel Jumpers with Felt Buttress
3134D	Coolcel Jumpers with Donut Buttress
3134W	Coolcel Jumpers - Wrap around (no buttress)



**Lateral J**

## 29 Knee Brace Undersleeve

This 2mm sleeve is to be worn under a more sophisticated knee brace to assist suspension. Also available in lightweight, breathable, latex-free Coolcel material.

**SIZING:** See chart on opposite page

Product #	Description
29	2mm Neoprene Undersleeve
29/3	3mm Neoprene Undersleeve
3129	Coolcel Undersleeve





### 3126 Infrapatellar Strap

**Indications:**

- Jumper's Knee (Patellar tendonitis), Osgood-Schlatter's

**Features:**

- Pain relief by increasing tension placed on the infrapatellar tendon - stabilizes the patella by reducing slack
- Rubber buttress concentrates the corrective force
- Easy to apply and adjust tension with Velcro closure
- Available in Coolcel in black

**SIZING:**

Size	Knee Center Circ.
S	11"-13"
M	13"-15"
L	15"-17"



### 6026 Patella Tendon Strap

**Features:**

- Consistent compression on the patella tendon to assist in correct patella tracking
- Reduces pain associated with a loose patella tendon
- Pressure maintained with a Griptech non-slip liner, holding the strap in place
- Adjustable straps ensure personalized fit and compression

**SIZING:** Universal



### 3511 Elastic Knee Support

**Indications:**

- General knee pain
- Arthritis

**Features:**

- Designed to provide warmth and mild support and compression
- Four-way breathable and seamless elastic design for easy knee movement

**SIZING:**

Size	Knee Circ.
XS	11"-14"
S	14"-15"
M	15"-16"
L	16"-17"
XL	17"-21"



### 3512 Charcoal Elastic Knee Support

**Features:**

- Designed to provide warmth, mild support, and compression
- Breathable, seamless elastic and nano bamboo material to inhibit bacteria and odour
- Four-way knit design allows for easy movement

**SIZING:**

Size	Knee Circ.
S	14"-15"
M	15"-16"
L	16"-17"
XL	17"-21"



### 5530 3D Elastic Knee Stabilizer

**Features:**

- 3D knitting technology provides superior fit and comfort
- Provides medical compression of 20-20mmHg to improve circulation and reduce swelling
- Spiral stays on both sides of the knee provide stability, while silicone dots on the inside stay covers stop the brace from sliding down

**SIZING:**

Size	Knee Circ.
S	11"-13"
M	13"-15"
L	15"-17"
XL	17"-19"
XXL	19"-21"

## 3611 Slip On Knee Compression

**Features:**

- Beige knitted elastic provides moderate compression
- 11" in length with knitted taper in the back

**SIZING:**

Size	Knee Circ.
S	12" - 14"
M	14.5" - 17"
L	17.5" - 20"
XL	20.5" - 24"



## 820 Sportec Knee Compression Sleeve

**Indications:**

- Strained knee joints
- Unstable muscle tissue

**Features:**

- Medical compression of 20-30mmHg
- Flexible spiral stays on the medial and lateral sides

Size	Knee Circumference
S	30-34cm
M	34-38cm
L	38-42cm
XL	42-46cm
XXL	46-50cm



## 3531 Sportec Patella Compression Knee Support

**Features:**

- Medical compression of 20-30mmHg
- Oval-shaped gel patella buttress stabilizes the patella and massages swelling away from the front of the knee
- Advanced 3D knitting ideal for active use, extremely comfortable for long-term wear
- Provides superior anatomical fit of the popliteal and muscle tissue for extreme support
- Flexible spiral stays on the medial and lateral sides provide extra stability

**SIZING:**

Size	Knee Circ.
S	30-35 cm
M	35-40 cm
L	40-45 cm
XL	45-50 cm



## 5537 Dynamic Knee Support

**Features:**

- 4 lateral stays for knee stability
- Gel patella pad for good patella support
- Cotton lining and breathable elastic material for comfort
- Contoured design for good fit and performance
- Velcro® compatible nylon provides universal size
- Wrap design for easy application

**SIZING:** Universal





### 3623 Spiral Elastic Knee Support

**Features:**

- Knitted beige elastic
- Tapered in the back for better fit
- Open patella with medial/lateral spiral stays for added stability

**SIZING:**

Size	Knee Circ.
S	11.5" - 13"
M	13" - 14.75"
L	14.75" - 16.5"
XL	16.5" - 17.25"



### Aspen ROM Knee

**Features:**

- Range of motion brace providing optimal compression and support
- Premium AirX material provides a high degree of breathability and moisture wicking, allowing for extended wear
- Also available in universal sizing upon request, made with lightweight neoprene material

**SIZING:** Universal (648111) or sized

Product #	Size	Knee Circ.
648101	XS	11 - 15"
648102	S	15 - 18"
648103	M	18 - 21"
648104	L	21 - 24"
648105	XL	24 - 27"
648106	XXL	27 - 35"



### Aspen OA Knee Wrap

**Features:**

- Decrease knee pain and inflammation associated with osteoarthritis
- Fits comfortably and securely to allow patients to return to active life
- Dual hinge feature ensures even pressure distribution
- Allows for flexion/extension stops

**SIZING:** 6" above and below knee centre


Product # Left	Product # Right	Size	Thigh Circ.	Calf Circ.
651021	651026	S	15.5 - 18.5"	12 - 14"
651022	651027	M	18.5 - 21"	14 - 16"
651023	651028	L	21 - 23.5"	16 - 18"
651024	651029	XL	23.5 - 26.5"	18 - 20"
651025	651030	XXL	26.5 - 29.5"	20 - 22"

### 150140 Bort Stabilo Knee Support with Hinge

Circular knitted elastic knee support for superior fit and comfort with side joint splints for guidance of the knee joint. Compression class II.

Features:

- The extension limit is 0 degrees
- Circular fastening straps enable the hinges to provide extra lateral stabilization to the knee joint
- Pressure reduction zones in the extremities and popliteal fossa areas improve function and fit
- Moisture wicking fabric and antimicrobial finish wicks perspiration away, keeping the knee cool during activity and eliminating odour
- The top and bottom edges are made of extra soft snag proof material for extra comfort



**SIZING:**

Size	Calf Circ. (d)	Thigh Circ. (f)
S	26 - 29 cm	32 - 36 cm
M	29 - 33 cm	36 - 40 cm
L	33 - 37 cm	40 - 44 cm
XL	37 - 41 cm	44 - 48 cm

Also available in plus sizes



### 114420 Bort Select Stabilogen Knee Support



Circular knitted elastic knee support for superior fit and comfort with silicone patella fixation support for stabilization and relief of knee joint pain. Compression class II.

Features:

- Two lateral stays on each side provide lateral support
- The DynamicFlex functional knitted fabric relax zone at the knee centre relieves and protects flexor tendons and soft tissues
- The lengthwise elastic stretch zone above the silicone pad ensures an ideal fit at all bending angles
- Moisture wicking fabric and antimicrobial finish wicks perspiration away, keeping the knee cool during activity and eliminating odour
- The top and bottom edges are made of extra soft snag proof material for extra comfort

**SIZING:** See 150140 Bort Stabilo Knee Support with Hinge



### 114700 Bort Asymmetric Knee Support



Circular knitted elastic knee support for superior fit and comfort with individual guidance of the patella in the medial and distal direction by means of a silicone half ring supported by lateral elastic straps. Compression class II.

Features:

- One lateral stay on each side provides lateral stability
- A silicone adhesive border provides a secure fit
- The DynamicFlex functional knitted fabric relax zone in the knee center relieves and protects flexor tendons and soft tissues
- Lengthwise elastic stretch zone in the lateral area of the silicone pad improves function and fit
- Moisture wicking fabric and antimicrobial finish wicks perspiration away, keeping the knee cool during activity and eliminating odor
- The top and bottom edges are made of extra soft snag proof material for extra comfort

**SIZING:** See 150140 Bort Stabilo Knee Support with Hinge





### 6640 Patella Stabilizer

**Features:**

- Two removable patella C-shaped buttresses to support the patella
- Removable, lightweight, flexible lateral stays improve knee stability
- Lycra popliteal area for maximum air circulation
- Lace elastic strap closure system for adjustability

**SIZING:**

Size	Thigh Circ. (10cm above knee centre)
Small	32 - 36cm
Medium	36 - 40cm
Large	40 - 45cm
XLarge	44 - 49cm



### 6642 Multi-Stability Lace Knee Brace

**Features:**

- Comes with polycentric hinge stays and spiral stays for variable levels of support
- Two elastic straps provide adjustable compression
- Soft silicone patella buttress for patella stability
- Stays provide lateral stability to support the MCL-LCL
- 100% cotton in contact with the skin for breathability and comfort
- Specify left or right

**SIZING:**

Size	Thigh Circ. (10cm above knee centre)
Small	32 - 36cm
Medium	36 - 40cm
Large	40 - 45cm
XLarge	44 - 49cm
XXL	48 - 54cm
XXXL	54 - 60cm
XXXXL	60 - 66cm



### 6622 Hinge Knee Wrap

**Features:**

- Silicone ring provides patella stability
- Removable, mouldable polycentric hinges provide medial and lateral support to MCL
- Lycra popliteal area provides circulation
- Opens completely for easy application with two elastic straps to adjust support/compression

**SIZING:**

Size	Thigh Circ. (10cm above knee centre)
Small	35 - 40cm
Regular	40 - 48cm
Extra	48 - 58cm



### 6653 ACL-PCL/MCL-LCL Stability Knee Brace

**Features:**

- Two shaped stays (medial and lateral) provide knee stability
- Variable rotation center to mimic natural knee movement, both extension and flexion can be limited to desired knee angle
- Easy to wear with adjustable straps and padding
- Special 3D Spacer woven fabric for breathability and comfort
- Specify left or right

**SIZING:**

Size	Thigh Circ. (10cm above knee centre)
XSmall	34 - 38cm
Small	38 - 43cm
Medium	43 - 49cm
Large	49 - 56cm
XLarge	56 - 64cm
XXL	64 - 73cm
XXXL	73 - 83cm

### 5438 Post-Op Knee Brace

**Indications:**

- Post surgery care of ACL, PCL, and other surgical repairs

**Features:**

- Trimmable foam wrap and adjustable straps allow for universal sizing and comfortable wear
- Easy to change dials can be set to a desired extension and flexion angle or set to a locked position
- Controls range of motion of the knee to protect newly constructed grafts, optimize ligament healing, enhance knee stability, and reducing the risk of re-tear
- Offers security and confidence to resume activity

**SIZING:**

Product #	Leg Length
5438S	16-19"
5438L	20-23"



### Universal 3-Panel Knee Immobilizer

**Indications:**

- Post knee surgery
- Following patellar dislocation
- Certain knee fractures

**Features:**

- Plush foam exterior
- Two medial and two lateral malleable stays
- Three posterior movable side stays and front Velcro closure to provide universal fit
- Popliteal foam pad for better fit and comfort
- Commonly referred to as a Zimmer splint

**SIZING:** Universal

Product #	Description
122	14" Length
123	20" Length
124	24" Length
125	22" Length



### 11-4 Glyde OA

**Features:**

- Simple, effective, and elegant approach to reducing osteoarthritis knee pain
- One strap, with an easy push button closure suspends the brace over the calf and provides dynamic pressure relief
- Breathable DiamoCool fabric wrap provides comfortable compression and conceals the 3-point leverage corrective frame

**SIZING:**

Size	6" above mid-patella	6" below mid-patella
S	14" - 17"	11.5" - 13.5"
M	17" - 20"	13.5" - 15.5"
L	20" - 24"	15.5" - 18"
XL	24" - 28"	18" - 20.5"
XXL	28" - 32"	20.5" - 23.5"



### 06-5 TRUalign OA

**Features:**

- Single upright OActive brace design in a patient-preferred wrap-style brace
- Q-hinge technology that provides angular correction to relieve the pressure and pain of knee osteoarthritis
- Easy to apply, breathable, comfortable, and conforms to anatomy
- Specified left and right

**SIZING:**

Size	6" above patella	6" below patella
S/M	14" - 18.5"	11" - 14"
M/L	18.5" - 25.5"	14" - 19"
L/XL	23.5" - 32"	19" - 25"





### Aspen OA Knee+

**Features:**

- Premium brace is ideal for patients with mild to moderate osteoarthritis
- Low profile, universal design allows for less inventory and better patient compliance
- Available for the left and right side

**SIZING:**

Product #	Side
651011	Left
651012	Right

**SIZING:** Universal



### Aspen CI Knee

**Features:**

- Dual strut ligament knee brace to address multiple types of instability
- Help improve confidence with increased mobility
- Flexion/extension stops can be added to the hinge to control allowed range of motion:
  - Extension stops range from 0 - 25 degrees
  - Flexion stops range from 0 (180) - 90 degrees
- Left or right specified

**SIZING:**

Product # Left	Product # Right	Size	Thigh Circ.	Calf Circ.
652001	652006	S	39-47 cm	30.5-35.5 cm
652002	652007	M	47-53.25 cm	35.5-40.5 cm
652003	652008	L	53.25-59.5 cm	40.5-47 cm
652004	652009	XL	59.5-67.25 cm	47-50.75 cm

For added comfort and to help reduce brace migration, try the suspension wrap strap, available with code **651111**.





## Sportec Compression Sleeves

Sportec Compression Sleeves are designed with medical grade compression, providing maximum pressure at the distal extremity and progressively decreasing to improve circulation and remove lactic acid from the muscles. Increased circulation delivers more oxygen to the muscles for maximum endurance and recovery. Specialized knitting avoids pressure points and chafing, providing superior comfort.

### 811 Sportec Calf Compression Sleeve

**Features:**

- Medical grade compression of 20mmHg; maximum pressure at the ankle, progressively less pressure to the top of the calf
- Fine elastic top band provides durable anti-slip with maximum comfort

**SIZING:** Ankle and calf circumference (in cm)



Size	cB	cC
XS	17-20	26-32
S	20-23	30-36
M	23-26	34-40
L	26-29	38-44
XL	29-32	42-48



### 805 Sportec Ankle/Calf Compression Sock

**Features:**

- Medical grade compression of 20mmHg; maximum pressure at the foot and ankle, progressively less pressure to the top of the calf
- Targeted zones support the achilles tendon and calf muscle to relieve muscle pain
- Special knitting supports arch, reinforced cushion sole absorbs impact and shock, loose-knit toe prevents blisters
- Nano-silver antimicrobial fibre provides odour reduction

**SIZING:** Ankle and calf circumference (Comes in pairs)



Size	Ankle Circumference	Calf Circumference
S	19 - 21.5cm	28 - 38cm
M	21.5 - 24cm	33 - 43.5cm
L	24 - 26.5cm	35.5 - 46cm
XL	26.5 - 29cm	38.5 - 51.5cm



### 815 Sportec Leg Compression Sleeve

**Features:**

- Special knitting technique that targets the calf and thigh muscles to enhance strength and energy
- Reinforced knitting surround the patella to provide mild support for injury prevention,
- Compresses tendons and ligaments to enhance strength in muscle contraction and stretch
- Thin breathable construction efficiently transports body heat and perspiration away from skin
- Comes in singles

**SIZING:** Calf and thigh circumference (in cm)



Size	XS	S	M	L	XL	XXL
cX	40-48	44-52	48-56	52-60	56-64	60-68
cY	25-29	29-33	33-37	37-41	41-46	46-51
cZ	17-20	20-23	23-26	26-29	29-32	32-35





Product No.	Product	Page
-------------	---------	------

### Ankle Supports - Swede-O

5172	Swede-O Ankle Lok	H2
5173	Swede-O Strap Lok	H2
5179	Swede-O Inner Lok 8	H2
5180	Swede-O Tarsal Lok	H2

### Ankle Supports

5570	Dynamic Ankle Lacer	H3
5370	Active Strap Lock Ankle Brace	H3
5375	Air Cushion Ankle Brace	H3
74B	Active Ankle Trainer II	H3
1450	Bort ActiveColour Ankle Support	H4
1400	Slip-On Ankle Compression	H4
325	Stromgren 2 Strap Ankle Support	H4
054100SP	Bort Sport Gel Ankle	H4
1404	Elastic Ankle Wrap	H5
5192	FXT Compression Socks	H5
801	Sportec Compression Socks	H5
805	Sportec Ankle/Calf Compression Sleeve	H5

### Ankle Sleeves

R2622	Dynagel Heel Protector	H6
35038/35039	Intelligel Performance Achilles Protector	H6
35084	Intelligel Performance Shin Guard	H6
35186	Intelligel Performance Boost Bumper	H6

Product No.	Product	Page
-------------	---------	------

### Walking Boots

5354	Pro Ankle Walking Boot	H7
5355	Pro Ankle Walking Boot Short	H7
5358	Ultralight Air Walking Boot	H7
5359	Ultralight Air Walking Boot Short	H7

### Ankle Foot Orthoses

71	Toe Lift	H7
1600	Swedish Ankle Foot Orthosis	H8
6970	Sky Medical Posterior Leaf Spring AFO	H8
6975	Sky Medical Semi-Solid AFO	H8
13510-13518	Orthomerica Leaf Spring Orthosis (LSO)	H9
13540-13548	Orthomerica Ankle Foot Orthosis (AFO)	H9
13500-13508	Orthomerica Supra Lite AFO	H9
13520-13528	Orthomerica Carbon Supra Lite AFO	H9
8000	Easy Walk Advance AFO	H10
8010	Easy Walk Eco AFO	H10
8012	Easy Walk Advance II AFO	H10
8020	Easy Walk Stability AFO	H11
8030	Easy Walk Fixation AFO	H11
8100	Step-On 100 Lateral Hinge AFO	H12
8101	Step-On 100 Medial Hinge AFO	H12
8110	Step-On 200 Dual Hinge AFO	H12
FH250	Trulife Matrix Max AFO	H13
FH260	Trulife Matrix Max 2 AFO	H13
FH270	Trulife Matrix SuperMax AFO	H13



## 5172 Ankle Lok

**Features:**

- Triple layer, woven vinyl-laminate support featuring “Ankle Lok” offset flap which tightly pulls the ankle into a stable position and secures it when laced
- U-shaped, internal stainless steel spiral stays provide extra support
- Removable plastic stabilizers can be slid into retention tabs to give maximum support

**SIZING:** Shoe size

Size	XS	S	M	L	XL	XXL
Men	3 - 5	6 - 7	8 - 10	11 - 12	13 - 14	15+
Women	4 - 6	7 - 9	10 - 11	12 - 13	--	--



## 5173B Strap Lok

**Features:**

- Lightweight, non-stretch nylon straps wrap in figure eight configuration, locking the ankle in protected position
- The strap at the top of the brace holds the figure eight in place and minimizes lace movement and loosening
- Available in black

**SIZING:** Shoe size

Size	Men	Women
XS	3 - 5	4 - 6
S	6 - 7	7 - 9
M	8 - 10	10 - 11
L	11 - 12	12 - 13
XL	13 - 14	--
XXL	15+	--



## 5180 Tarsal Lok

**Features:**

- Combines the support of a rigid ankle brace with the comfort of a lace-up brace
- Built-in stabilizers mould to the shape of the ankle from body heat
- Low profile design allows it to be worn in almost any shoe and with orthotics
- Speed lace system allows for quick and easy brace application
- Available in black

**SIZING:** Shoe size

Size	Men	Women
XS	3 - 5	4 - 6
S	6 - 7	7 - 9
M	8 - 10	10 - 11
L	11 - 12	12 - 13
XL	14+	--



## 5179 Inner Lok 8

**Features:**

- Two figure eight straps are pre-configured through the inside of the brace in the correct position to ensure proper strap placement every time
- Straps provide extra support and allow for a more custom fit
- Patients can apply the brace with very little instruction with quick and easy application

**SIZING:** Shoe size

Size	S	M	L	XL	Product #	Description
Men	6 - 7	8 - 10	11 - 12	14+	<b>5179</b>	White Inner Lok 8
Women	7 - 9	10 - 11	12 - 13	--	<b>5179B</b>	Black Inner Lok 8

### 5570 Dynamic Ankle Lacer

**Features:**

- Laces with Velcro tabs for quick and easy application and tightening
- Incorporates breathable fabric for cool comfort and bilateral plastic stays for added ankle stability

**SIZING:**

Size	Ankle Circ.
S	16 - 21 cm
M	21 - 26 cm
L	26 - 31 cm
XL	31 - 37 cm



### 5370 Active Strap Lock Ankle Brace

Designed to stabilize the active ankle in sports activities to prevent injury or to stabilize a previously injured ankle.

**Features:**

- Quick-lace system allows for easy application and removal
- Plastic side bracing system provides optimal support to limit ankle inversion and eversion while allowing flexion and extension for sport activities
- Figure eight strapping system provides added stability needed for high-activity use

**SIZING:**

Size	Ankle Circ.
XS	8-9"
S	9-10"
M	10-11"
L	11-12"
XL	12-13"
XXL	13-14"



### 5375 Air Cushion Ankle Brace

**Indications:**

- Acute ankle injury, chronic ankle instability/weakness

**Features:**

- Supports the ankle after sprain for pain reduction
- Heel tab on the bottom of the brace can be adjusted for the width of the foot
- Velcro can be tightened to the correct size for the ankle
- Air cushions come pre-inflated for immediate use and a small pump is included to re-inflate the cushions if needed

**SIZING:** Universal



### 74B Active Ankle Trainer II



**Features:**

- "T2" features a durable single strap system that is adjustable for different shoe heights
- Feather-light EVA padding hugs the contours of the ankle
- Sleeker, less bulky brace while still offering maximum performance and protection
- Provides medial/lateral support and compression to the joint

**Indications:**

- 1st or 2nd degree ankle sprains
- Chronic or recurrent problems as well as instabilities in the ankle
- Ideal for ankle protection while playing basketball, football, or volleyball

**SIZING:** Shoe size

Size	Men	Women
S	5 - 8.5	6.5 - 9.5
M	9 - 12	10 - 13
L	12.5 - 16	13.5 - 16





## 1450 Bort ActiveColor Ankle Support



**This item has been discontinued, limited stock available. Please contact to confirm.**

Provides a perfect fit with improved edge and reciprocated heel. Offers effective bracing and can be worn with sport and walking shoes. Comes in black.

**Indications:**

- Overstrain prophylaxis
- Sprains and strains
- Irritations caused by overstraining

**SIZING:**

Size	Ankle Circ.
S	up to 21 cm
M	21-23 cm
L	23-25 cm
XL	25-27 cm
XXL	27-29 cm



## 1400 Slip On Ankle Compression

**Features:**

- Beige, knitted elastic provides intermediate compression
- Open toe and heel
- Knitted taper for contoured fit

**SIZING:**

Size	Ankle Circumference
S	7" - 8"
M	8.25" - 9"
L	9.25" - 10"
XL	10.25" - 11.25"



## 325 Stromgren Double Strap Ankle Support

**This item has been discontinued, limited stock available. Please contact to confirm.**

**Features:**

- Beige Spandex that provides ankle/arch support
- Eliminates bulk and reduces tape use
- Double strap offers compression to medial and lateral sides of the ankle
- Velcro strap gives uniform compression/tension and allows breathability
- Wear over a sock, completely washable

**SIZING:** Specify left or right

Size	Ankle Circumference
XS	6" - 7"
S	7.5" - 8"
M	8" - 9.5"
L	9.5" - 11"
XL	11" - 12.5"



## 054100SP BORT Sport Gel Ankle

**This item has been discontinued, limited stock available. Please contact to confirm.**

**Features:**

- Compression of soft tissues, with Stabilo Strap, elastic can be freely positioned
- SOFTflex 3D knitted fabric with a velvety soft surface for more elasticity, stability, and comfortable wear
- Two silicone pads for intermittent passage
- COOLMAX wicks moisture away from the body for dryness and comfort
- Latex-free fabric

**SIZING:** Specify left or right

Size	Ankle Circumference
XS	up to 19cm
S	19 - 21cm
M	21 - 23cm
L	23 - 25cm
XL	25 - 27cm

### 1404 Elastic Ankle Wrap

**This item has been discontinued, limited stock available. Please contact to confirm.**

**Features:**

- 3" wide comfortable elastic
- Elastic wraps in a figure-eight style
- Hook and loop closure

**SIZING:**

Size	Ankle Circumference
S	6.5" - 7.5"
M	7.5" - 8.5"
L	8.5" - 10"
XL	10" - 12"



### 5192 FXT Compression Socks



**Features:**

- Targeted compression zones to aid in the pain management and relief from plantar fasciitis, heel spurs, and arch pain
- Ideal for sports or those on their feet all day

**SIZING:**

Shoe Size	XS	S	M	L	XL
Men's	4 - 5.5	6 - 7.5	8 - 9.5	10 - 11.5	12 - 14
Women's	5 - 6.5	7 - 8.5	9 - 10.5	11 - 12.5	13 - 15



### 801 Sportec Compression Socks

**Features:**

- Maximum compression level of 30mmHg
- Targeted achilles and plantar fasciitis compressions zones provide extra support of the achilles and foot arch
- Nano-silver antimicrobial fibre provides odour reduction
- Loose-knit toe prevents blisters
- Moisture-wicking construction for cool performance

**SIZING:**

Size	Ankle Circumference
XS	17 - 19cm
S	19 - 21.5cm
M	21.5 - 24cm
L	24 - 26.5cm
XL	26.5 - 29cm



### 805 Sportec Ankle/Calf Compression Sock

**Features:**

- Medical grade compression of 20mmHg; maximum pressure at the foot and ankle, progressively less pressure to the top of the calf
- Targeted zones support the achilles tendon and calf muscle to relieve muscle pain
- Special knitting supports arch, reinforced cushion sole absorbs impact and shock, loose-knit toe prevents blisters
- Nano-silver antimicrobial fibre provides odour reduction

**SIZING:** Ankle and calf circumference (Comes in pairs)



Size	Ankle Circumference	Calf Circumference
S	19 - 21.5cm	28 - 38cm
M	21.5 - 24cm	33 - 43.5cm
L	24 - 26.5cm	35.5 - 46cm
XL	26.5 - 29cm	38.5 - 51.5cm





**A**



**C**



**B**



**D**

**(A) R2622 Dynagel Heel Protector**

**This item has been discontinued, limited stock available. Please contact to confirm.**

**Features:**

- Comfortable gel pad in the back of sleeve
- Reduces pressure and friction where the Achilles attaches to the heel
- Anatomically-shaped open-toe elastic sleeve designed to hold the gel pad in the correct position
- One per box

SIZING: Universal

**(B) Intelligel Performance Achilles Protector**

**This item has been discontinued, limited stock available. Please contact to confirm.**

**Features:**

- Lightweight, comfortable toeless sleeve
- With no bulk, can be worn with casual and dress shoes
- Anatomically-shaped, soft gel pad comforts, cushions, and protects the tender achilles area

**Indications:**

- Achilles tendonitis and bursitis
- Heel prominences
- Hardened or cracked skin

Product #	Size
35038	S/M
35039	L/XL

**(C) 35084 Intelligel Performance Shin Guard**

**This item has been discontinued, limited stock available. Please contact to confirm.**

**Features:**

- 3" by 5" elastic tube
- Lined with gel to reduce pressure, pinching, or friction on the lower leg
- Ideal for use in ice skates or ski boots
- One per box

SIZING: Universal

**(D) Intelligel Performance Boot Bumper**

**This item has been discontinued, limited stock available. Please contact to confirm.**

**Features:**

- Conforming 1/8" polymer gel pads protect the front and back of the ankle from pressure, pinching, and friction
- Intelligel Elastic sleeve soothes irritation, burning, and redness
- Protects the ankle from abrasion caused by athletic footwear without added bulk - ideal for use in ice skates or ski boots

Product #	Size
35186	S
35187	L



### 5354/5355 Pro Ankle Walking Boot

This low-profile foam liner walking boot can be used instead of a cast for the treatment of stable fractures and ankle sprains. The lightweight rocker sole design reduces shock transfer during heel strike, and is engineered to promote a natural, stable gait. The contoured struts are designed to provide maximum patient comfort. The foot section can accommodate full dressing and padding. Maximizes strength, reduces overall weight, and minimizes bulkiness. Fits left or right. Available in short and long.

**Indications:**

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable stress fractures
- Post-operative use
- Trauma and rehab

**5354 & 5355 SIZING:** Shoe size

Size	Men's	Women's
XS	< 4	< 5.5
S	4-7	6-8
M	7-10	8-11
L	10-12	11-13
XL	12+	13+

Product #	Description
5354	Walking Boot Regular
5355	Walking Boot Short



### 5358/5359 Ultralight Air Walking Boot

The Ultralight Air Walking Boot has an inflatable foam liner that allows for simple, tool-free adjustment to achieve custom, fine-tuned, consistent compression. Straps fasten in the front for easy application, adjustment and removal. The lightweight sole and cushioned heel provide excellent shock absorption, and the rocker bottom is low and wide to help promote a natural, stable gait. Ventilated panels reduce microbial buildup.

**Indications:**

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable fractures
- Post-operative use
- Trauma and rehab

**5358 & 5359 SIZING:** Shoe size

Size	Men's	Women's
XS	< 4	< 5.5
S	4-7	6-8
M	7-10	8-11
L	10-12	11-13
XL	12+	13+

Product #	Description
5359	Ultralight Air Walking Boot Short
5358	Ultralight Air Walking Boot Regular



### 71 Toe Lift

**This item has been discontinued, limited stock available. Please contact to confirm.**

The toe lift helps provide toe clearance during walking. A low profile plastic attachment as well as a metal attachment for shoes with narrow eyelet spacing are provided with this support. A slipper attachment is also provided for wearing in the house without shoes.

**Indications:**

- Drop foot
- Toe dragging while walking

**SIZING:** Shoe size, ankle circumference

Size	Women's	Men's	Ankle Circumference
S	5-8	6	up to 10"
M	9-11	7-10	up to 11"
L	12-14	10-14	up to 14"
XL	-	15-17	14-18"





## 1600 Swedish Ankle Foot Orthosis

Effective for patients suffering from a variety of foot ailments, providing static dorsiflexion assistance and lateral stability for the entire foot/ankle area. Quick, efficient trimming to size with scissors, eliminating the need for a heat gun. The low arch and open heel give a streamlined profile that fits easily into any shoe. Injection moulded polyethylene allows for different thicknesses throughout the orthosis - thicker on the vertical aspect for rigidity and a thinner footplate for trimmability. Padding is provided on the calf section for extra comfort and a Velcro® ankle strap secures the Swedish AFO around the ankle. The orthosis' height is 11".

**SIZING:** Universal

Product	Description	Product	Description
1600LR	Left Regular	1600LX	Left XL
1600RR	Right Regular	1600RX	Right XL

## 6970 Posterior Leaf Spring AFO



This lightweight AFO is designed to assist ambulation for those patients suffering from mild to moderate neurological deficit, resulting in drop-foot. The anatomically correct shape will fit inside most lace-up shoes. It is fabricated from specially-formed, orthopedic-grade thermoplastic, and offers durable, long-term performance and reliability. For a more comfortable individual fit, the orthosis may be custom trimmed and reshaped using cast scissors and a heat gun.

**SIZING:** Shoe size, specify left or right

Size	S	M	L	XL
Men	5 - 6	6.5 - 10	10.5 - 13	13 +
Women	6.5-7.5	8-11.5	12-14.5 +	14.5 +

## 6975 Semi-Solid AFO



This robust AFO is designed to provide increased dorsiflexion assist and greater medial-lateral control during ambulation for those patients suffering neurological deficit resulting in drop-foot. The anatomical shape will fit inside most lace-up shoes. It is fabricated from specially-formed orthopedic grade thermoplastic and offers durable long term performance and reliability. For a more comfortable individual fit, the orthosis may be custom trimmed and reshaped using cast scissors and a heat gun.

**SIZING:** Shoe size, specify left or right

Size	S	M	L
Men	5 - 6	6.5 - 10	10.5 - 13
Women	6.5-7.5	8-11.5	12-14.5 +

## Orthomerica Leaf Spring Orthosis

Popular design accommodates patients not requiring a custom or more rigid AFO. For flaccid foot drop indications. Fits in most normal lace-up shoes. Low profile for excellent cosmesis. Can be easily trimmed and heated for a more custom fit.

**SIZING:** Shoe size, specify left or right

Size	Men	Women	Left	Right	Height
S	5 - 7	6 or less	13510	13515	30cm
M	7.5 - 10	6.5 +	13511	13516	34cm
L	10.5 - 13	-	13512	13517	34cm
XL	13+	-	13513	13518	36cm



## Orthomerica Ankle Foot Orthosis

Prefabricated from 1/8" polypropylene for necessary rigidity and comfort. Indicated for post stroke flaccid foot drop or drop foot spasticity. Can be heat treated for further modification and individual fit. Very versatile design.

**SIZING:** Shoe size, specify left or right

Size	Men	Women	Left	Right	Height
S	5 - 7	6 or less	13540	13545	30cm
M	7.5 - 10	6.5 +	13541	13546	34cm
L	10.5 - 13	-	13542	13547	34cm
XL	13+	-	13543	13548	36cm



## Orthomerica Supra Lite AFO

Comfortable, low profile prefabricated leaf spring design to correct flaccid foot drop, especially on post-stroke patients. Injection moulded polypropylene allows for a thicker on the vertical aspect for rigidity and thinner footplate for ease of trimming, comfort and flexibility. Special heel cutout accommodates a wide range of men's and women's sizes. \*Maximum weight of 180 lbs.

**SIZING:** Shoe size, specify left or right

Size	Men	Women	Left	Right	Height
S	5 - 7	6 or less	13500	13505	10.5"
M	7.5 - 10	6.5 +	13501	13506	11.5"
L	10.5 - 13	-	13502	13507	13.5"
XL	13 +	-	13503	13508	15.5"



## Orthomerica Carbon Supra Lite AFO

This injection-moulded, carbon-reinforced polypropylene AFO designed to fit patients requiring an AFO for flaccid foot drop. Lightweight, comfortable design with carbon rebound feature provides dorsiflexion, and the open heel allows it to fit in most shoes. Padded tibial crest and dorsum straps enable accurate and secure positioning. Please note, this is not a solid carbon AFO.

**SIZING:** Shoe size, specify left or right

Size	Men	Women	Left	Right	Height
S	5 - 7	6 or less	13520	13525	10.5"
M	7.5 - 10	6.5 +	13521	13526	11.5"
L	10.5 - 13	-	13522	13527	13.5"
XL	13 +	-	13523	13528	15.5"





## 8000 Easy Walk® Advance AFO

Made of pure carbon fibre. The construction of the orthosis in the calf area allows limited ankle joint movement. The sole part near the toes and the heel is more flexible than in the middle. This together with usual for the line “spring” property make walking with the orthosis less tiring. Due to an especially thin sole, the orthosis does not affect the posture.

This AFO is prefabricated and requires individual adjustment. It can be fitted by either grinding or simply cutting the sole with scissors.

**SIZING:** Shoe size, specify left or right

Size	Men's	Women's	Height
XS	--	4 - 6	33 cm
S	5 - 7	6 - 9	35 cm
M	7 - 10	9 - 12	36 cm
L	10 - 13	12+	37 cm
XL	13 - 16	--	39 cm



## 8010 Easy Walk® Eco AFO

Produced by combing carbon fibre with composite materials. Comparison tests have confirmed that the product kept the functionality and comfort level close to pure carbon fibre version, with the exception of slight difference in weight. Due to lower resilience of composite materials, the “spring” effect here is slightly softer if compared with pure carbon orthoses and that could serve as advantage for some patients.

This AFO is prefabricated and requires individual adjustment. It can be fitted by either grinding or simply cutting the sole with scissors.

**SIZING:** Shoe size, specify left or right

Size	Men's	Women's	Height
XS	--	4 - 6	33 cm
S	5 - 7	6 - 9	35 cm
M	7 - 10	9 - 12	36 cm
L	10 - 13	12+	37 cm
XL	13 - 16	--	39 cm



## 8012 Easy Walk® Advance II AFO

Made from pre-preg carbon and optimized for the toe-lift while staying lightweight and thin. Calf support strut is designed to adapt to the variety of shapes and provide energy return during the gait. New strapping system with convenient magnetic clasp combines firm hold, durability and easily adjustable length.

This AFO is prefabricated and requires individual adjustment. It can be fitted by either grinding or simply cutting the sole with scissors.

**SIZING:** Shoe size, specify left or right

Size	Men's	Women's	Height
XS	--	4 - 6	32.5 cm
S	5 - 7	6 - 9	33.5 cm
M	7 - 10	9 - 12	34.5 cm
L	10 - 13	12+	35.5 cm
XL	13 - 16	--	37.5 cm

### 8020 Easy Walk® Stability AFO

Made by combing carbon fibre with composite materials. It is designed to provide dynamic assistance with ankle instability control. With the three-point fixation system (at the patient's calf, ankle and foot), medial-lateral as well as rotational control of the foot and ankle complex is fully handled. The specific design of the support strut in the "stability" AFO allows maximum force distribution and energy return during the gait.

This AFO is prefabricated and requires individual adjustment. It can be fitted by either grinding or simply cutting the sole with scissors.

**SIZING:** Shoe size, specify left or right

Size	Men's	Women's	Height
XS	--	4 - 6	32.5 cm
S	5 - 7	6 - 9	33.5 cm
M	7 - 10	9 - 12	34.5 cm
L	10 - 13	12+	35.5 cm
XL	13 - 16	--	37.5 cm



### 8030 Easy Walk® Fixation AFO

Made from pure carbon fibre is the newest model in the line. It is designed to maximize ankle support and provide optimal biomechanical gait cycle. This dynamic ankle-foot orthosis allows patient to have medial-lateral control of the foot and ankle complex due to strong lateral fixation of lower leg. Due to the specific design of the support strut that allows maximum lateral support during the gait, major usual for pre-preg items energy return in this orthosis is carried out by the flexible footplate.

This AFO is prefabricated and requires individual adjustment. It can be fitted by either grinding or simply cutting the sole with scissors.

**SIZING:** Shoe size, specify left or right

Size	Men's	Women's	Height
XS	--	4 - 6	32.5 cm
S	5 - 7	6 - 9	33.5 cm
M	7 - 10	9 - 12	34.5 cm
L	10 - 13	12+	35.5 cm
XL	13 - 16	--	37.5 cm



## Step-On AFO

Step-On AFOs are extremely lightweight and incorporate unilateral and bilateral spring-loaded hinges. This product line is unique because it offers off-the-shelf AFOs, yet provides articulation and adjustable dorsiflexion assistance which can be easily adjusted with an allen key to change the spring tension within the hinge. This functionality provides patients with the exact amount of control they need through every stage of rehabilitation, while still allowing essential functional movements during the gait cycle.



### 8100 100 Lateral Hinge AFO

**STEP-ON**

This brace is the go-to for most cases of drop-foot, a unilateral hinge located on the lateral side of the ankle for an easy fit into normal footwear. This AFO only weighs 290g.

**SIZING:** Specify left or right

Size	Shoe Size (Men)	Shoe Size (Women)
S/M	3-7.5	5-9.5
L/XL	8-12	10-13



### 8101 100 Medial Hinge AFO

**STEP-ON**

This brace is ideal for patients with uncorrectable varus and prominent lateral malleolus. It includes a unilateral hinge on the medial side of the ankle, and can easily fit in normal footwear. This AFO only weighs 290g.

**SIZING:** Specify left or right

Size	Shoe Size (Men)	Shoe Size (Women)
S/M	3-7.5	5-9.5
L/XL	8-12	10-13



### 8110 200 Dual Hinge AFO

**STEP-ON**

This brace is for cases of drop-foot with excessive correctable varus/valgus (supination/pronation) of the foot, featuring bilateral hinges and a patented strap system to provide effective, comfortable support. This AFO only weighs 458g.

**SIZING:** Specify left or right

Size	Shoe Size (Men)	Shoe Size (Women)
S	3-5.5	5-7.5
M	6-9	8-10
L	9.5-10.5	11.5-12.5
XL	11-13	13+

## Trulife Matrix AFOs

The Matrix line of AFOs from Trulife are made with an energy-storing composite footplate, an open calcaneal design, and a lateral strut. The foot plate can be trimmed for best fit and the padded anterior shell is height-adjustable. The carbon is warranted for one year and the softgoods for 90 days. Two tibial pads are included.

### FH250 Matrix Max



The Matrix Max is Trulife's best-selling Matrix AFO. The reinforced lateral strut, made of lightweight carbon composite, provides enhanced stability. Suitable for low to high activity patients with unilateral or bilateral drop foot, as well as partial foot amputations.

**SIZING:** Specify left or right

Product # - Left	Product # - Right	Size	Height	Foot Length
FH250L-02	FH250R-02	S	39-42cm / 15.5-16.5in	<25cm / <9.8in
FH250L-03	FH250R-03	M	41-43cm / 16-17in	24-27cm / 9.4-10.6in
FH250L-04	FH250R-04	L	43-46cm / 17-18in	25-29cm / 9.8-11.4in
FH250L-05	FH250R-05	XL	43-46cm / 17-18in	26-31cm / 10.2-12.2in



### FH260 Matrix Max 2



This AFO has all the same features as the Matrix Max. It also has an easy-trim foot plate that allows fitting of up to 3 shoe sizes with one product. It has a heat-mouldable, trimmable and height adjustable anterior shell. A narrower foot plate in the medial arch allows for an easier shoe fit. Shrink tubing along the lateral strut protects the shoe against damage.

**SIZING:** Specify left or right

Product # - Left	Product # - Right	Size	Height	Foot Length
FH260L-02	FH260R-02	S	39-42cm / 15.5-16.5in	<25cm / <9.8in
FH260L-03	FH260R-03	M	41-43cm / 16-17in	24-27cm / 9.4-10.6in
FH260L-04	FH260R-04	L	43-46cm / 17-18in	25-29cm / 9.8-11.4in
FH260L-05	FH260R-05	XL	43-46cm / 17-18in	26-31cm / 10.2-12.2in



### FH270 Matrix SuperMax



The Matrix SuperMax is the ultimate AFO, with a uniquely-designed footplate to provide greater energy retention at toe off. The shank angle is increased to 82° to improve knee extension movement at mid stance. Suitable for low to high activity patients with unilateral or bilateral drop foot, peroneal nerve palsy, incomplete paralyses, CVA or CMT.

**SIZING:** Specify left or right

Product # - Left	Product # - Right	Size	Height	Foot Length
FH270L-02	FH270R-02	S	39-42cm / 15.5-16.5in	<25cm / <9.8in
FH270L-03	FH270R-03	M	41-43cm / 16-17in	24-27cm / 9.4-10.6in
FH270L-04	FH270R-04	L	43-46cm / 17-18in	25-29cm / 9.8-11.4in
FH270L-05	FH270R-05	XL	43-46cm / 17-18in	26-31cm / 10.2-12.2in





Product No.	Product	Page
-------------	---------	------

**Foot Orthotics (OTS)**

80	Motion Control Arch Support	13
81	Novaped Insole	13
82	Semi-Rigid Arch Support	13
84	Novaped Full Length Rigid Insole	13
89	Active PF Foot Support	14
200-215	Novatherm Orthotics	14
224	Novaped Full Length Multiform Insole	14
213	Morton Extension Insole	15
26059702	Novaped Soft Orthotic Insoles	15
219/220	Novatherm Multi Heel Cup	15
295	Intercept Insole (Raised Heel)	15
204	Children's Orthotic - 3/4 length	16
205	Children's Orthotic - Full length	16
218	Children's Orthotic - Narrow	16
10761	Multiflex Reinforced Insole	16
10898	Plastoform Diabetic Insole	16

**Heel/Shock Absorbers**

HL	Adjust-A-Heel Lift	16
303	Intercept Heel Cushion	17
304	Intercept Heel Spur Pad	17
304C	Intercept Heel Spur Pad w/ Curl	17
R2673	Dynagel Comfort Heel Cup	17

**Formthotics**

306	Original Sport	18
307	Sport Cycle	18
308/309	Sport Ski/Ski Touring	18
310	Sport Football	18
311	Sport Run	18
312	Everyday Industry Workfit	18
313	Sport Hike Dual	18
316	ShockStop Diabetic Insole	18
314	Original Dual Medium	19
317	Original Single Layer	19
326	Low Volume Single	19
327	Low Profile Dual	19
339	Original Dual Hard	19
341	Comfort Single, 3/4 Length	19
328-331	3D Additions	110
318-324	Formthotics Orthotic Components	110
338	Wedge Strip	111
FHMV20	Formthotics Heater and Boot Extension	111

Product No.	Product	Page
-------------	---------	------

**Tuli's**

346/347	Tuli's Gel Heel Cup	112
351/358	Tuli's Classic/HD Heel Cup	112
352	So Soft Heel Cup	112
359	Tuli's Cheetah	113
357	Tuli's Metatarsal Cushion	113
349	Plantar Fasciitis Insole	113

**Miscellaneous Forefoot Supports**

R5580	Hallux Valgus Night Splint	114
479	Leather Metatarsal Pad - Kidney	114
20770	Leather Metatarsal Pad - Teardrop	114
112070	Metatarsal Bandage	114
810	Hallux Valgus Daytime Splint	114
26140	Hallux Valgus Night Splint	114
823	Toe Crest Cushion	114
829	Buddy Toe Splint Cushion	114
180	Hallux Valgus Daytime Splint	114

**Silicone Footcare**

CC206	Toe Separators	115
CC229	Elastic Band with Silicone Metatarsal Insert	115
CC220	Metatarsal Pad with Reinforced Loop	115
CC215	Heel Spur Cushion with Removable Cavity	115
CC218-2	Toe Spreader	115
CC256	Forefoot Cushions	115
CC226	Toe Crest with Reinforced Loop	115
CC227-2	Toe Caps	115
CC212	Heel Spur Cup	115

**Footcare**

20781	Pronation/Supination Wedge	116
20782	Heel Wedge	116
22400-22420	Longitudinal Arches	116
22500-22520	Metatarsal Bars	116
22600-22620	Metatarsal Pads	116
31150-31155	Toe/Finger Spreaders	116
31160-31163	Toe Separators	116
35038/35039	Intelligel Performance Achilles Protector	117
35186/35187	Intelligel Performance Boot Bumper	117
31030/31031	M-Gel Body Disk	117

Product No.	Product	Page
-------------	---------	------

### Dynagel Finger and Toe Gel Products

2600	Digital Tube	118
2601	Digital Cap	118
2604	Toe Separator	119
2605	Toe Spreader	119
2607	Bunion Sleeve	119
2608	Digital Pad	119
2609	Corn Pad	119
2615	Toe Spreader with Sleeve	119
2617	Toe Spreader/Bunion Protector	120
2622	Heel/Elbow Protector	120
2626	Digit Tube	118
2628	Toe Spreader with Two Loops	119
2629	Bunion Guard with Toe Spreader	120
2630	Bunion Guard	119
2631	Five Toe Separator	119
2632	Toe Crest	120
2633	Digital Cap	118
2634	5 <sup>TH</sup> Digit Valgus Guard	120
2636	Metatarsal Loop with Pad	120
2637	Toe Spreader with Loop	118
2647	Hammer Toe Splint	118

### Slippers and Compression Socks

801	Sportec Compression Socks	121
5192	FXT Compression Socks	121

### Plantar Fasciitis Night Splints

70U	Soft Plantar Fasciitis Night Splint	121
6070	PF Night Sock	122
73	Sole Stretch	122
5380	Plantar Fasciitis Adjustable Night Splint	122
6576U	Dorsal Night Splint	122

### Dynamic Offloading Shoes

595	Dynamic Post-Op Shoe	123
596	Dynamic Heel Offloading Shoe	123
597	Dynamic Forefoot Offloading Shoe	123

### Sizing Chart

Ortho Active	#1	#2	#3
<b>Metric</b>	35/36	37/38	39/40
<b>Women's</b>	5 - 6	7 - 8	9 - 10
<b>Men's</b>	--	5 - 6	7 - 8

Ortho Active	#4	#5	#6	#7
<b>Metric</b>	41/42	43/44	45/46	47/48
<b>Women's</b>	11 +	--	--	--
<b>Men's</b>	9 - 10	11 - 12	13 - 14	15 +

## 80 Motion Control Arch Support

**Features:**

- Support for the medial and lateral longitudinal arches and a cup for the heel
- Controls excess movement that causes foot pain
- Breathable top cover is skin-friendly surface that wicks perspiration
- 3/4 length enables fit into most shoes, can be used over glued-in insole
- Arch cavity on bottom allows fit over small arch support already in shoes
- Flexible EVA material is ideal for active use

**SIZING:** Shoe Size

Size	Men's	Women's	Size	Men's	Women's
1	--	5 - 6	4	9 - 10	--
2	5 - 6	7 - 8	5	11 - 12	--
3	7 - 8	9 - 10			



## 81 Novaped Insole

This ¾ length, medium density insole will fit easily into any shoe and provide support for both the metatarsal and longitudinal arches of the foot. 40 Shore A density Multiforte ensures good support with some flexibility for comfort and padding. The heel cup holds the foot in place while the flat bottom stops the insole from lateral rocking in the shoe.

Sizes 36-43

**SIZING:** See size chart on page J2



## 82 Semi-Rigid Arch Support

Semi-rigid orthotic with a rigid plastic longitudinal arch support, flexible heel and EVA top cover. Available with or without a metatarsal pad.

**SIZING:** Shoe size, see sizing chart for 80/83 above

Product #	Description
82	Regular
82M	With met pad



## 84 Novaped Full Length Rigid Insole

This full-length insole is made with the popular, ultra-thin Schein sandwich construction to provide a very thin and rigid support. It is fully finished so it can be worn just as it comes. It is also remoldable to your patient's cast if desired.

Sizes 36-43

**SIZING:** See size chart on page J2





## 89 Active PF Foot Support

The Active PF Foot Support is designed to support the active foot during sporting activities. The longitudinal arch support provides some flexibility, mimicking the natural motion of the arch to absorb shock during walking and running. The heel cup holds the foot in a stable position to eliminate unwanted foot motion. These support systems help control the movements that cause most types of foot pain. The padded top cover provides shock absorption during impact and the double padding in the forefoot provides extra comfort to the forefoot while running.  
**SIZING:** Sizes 1-6 (See size chart on page I-3)



200

## Novatherm Orthotic Insoles

Ortho Active brings you a versatile new line of German-designed prefabricated insoles. The Novatherm technology produces a very supportive, semi-rigid orthotic with a tremendously low profile. Easily customized with a heat gun, and can even be moulded to a cast. High flange and heel cup available for maximum control. Narrow cut version perfect for dress shoes. Order with topcover or apply your own, adheres well. Easily add postings or other modifications. Please see custom orthotic section for Novatherm orthotics without top cover.

**SIZING:** Comes in sizes 35-48 (see sizing chart on page I-3)



214

Product #	Sizes	Description
200	35-48	3/4 length with top cover, no met pad
206	35-48	3/4 length with top cover, met support
207	35-48	Full length with top cover, met support, high lateral wall and heel cup
214	36-40	3/4 length with top cover, met support and narrow fit
215	35-48	Full length with top cover, no met pad



215



## 224 Novaped Full Length Multiform Insole

Using their popular sandwich construction, Schein has created an insole that is very thin and includes a multiform top cover for long-lasting cushion comfort. The heel cup holds the foot in position, a solid arch provides excellent support, and a lateral flare holds the foot on top of the insole. Remould this insole to your patient's cast if desired for a custom fit and eliminate all the time spent cutting, gluing and finishing.

Sizes 36-43

**SIZING:** See size chart on page J2

## 213 Morton Extension Insole

This insole is a newer version of the **212**, a rigid insole ideal for patients with hallux rigidus, with a rigid base extending all the way to the first toe. This version features a comfortable top cover. Easily moulds to a patient's cast if desired.

**SIZING:** Sizes 36 - 45 (see size chart on page I-3)



## Novaped Soft Orthotic Insoles

The Novaped Soft Orthotic line is the perfect choice for a versatile, off-the-shelf insole, featuring a wide variety of products variations to suit the needs of your patients. Using Schein's signature sandwich construction, this line of orthotics is a quick solution that can be easily moulded to a cast to custom fit patients in a small amount of time, removing the steps of cutting, gluing and finishing.

Sizes 35-46

**SIZING:** See size chart on page J2

Product #	Top Cover	Sizes	Description
<b>26059702 - 26059713</b>	Black Fashmo with 25 Shore EVA	35 - 46	Original Novaped insole with Bircork for rigidity



## 219/220 Novatherm Multi Heel Cup

A full length prefab orthotic that can be remoulded. High heel counters are ideal for providing calcaneal stability. Consists of a pink or black multiform top, blue thermaforte bottom, and reinforcing middle layer.

**SIZING:**

Size	Children's Shoe Size	Size	Children's Shoe Size	Product #	Description
A	4 - 5	E	12 - 13	<b>219</b>	Multi Heel Cup Pink
B	6 - 7	F	1 - 2	<b>220</b>	Multi Heel Cup Black
C	8 - 9	G	3 - 4		
D	10 - 11	H	5 - 6		



## 295 Intercept Insole

Retail

Lightweight, high tech, polymer gel insole protects the foot from shock, absorbing and dissipating forces before they reach the body. Thin insole with a slightly raised heel, comes in sizes A-E.

**SIZING:** Shoe size

Size	Men	Women
A	2 - 3	4 - 5
A	3.5 - 4.5	5.5 - 6.5
B	5 - 6	7 - 8
B	6.5 - 7.5	8.5 - 9.5
C	8 - 9	10 - 11
C	9.5 - 10.5	11.5 - 12
D	11 - 12	
D	12.5 - 13.5	
E	14 - 15	
F	15.5 - 16.5	





205

## Suprinara Children's Foot Orthotics

Suprinara construction with standard lateral support and deep heel cup. Heat will not transfer through multiform topcover, allowing the clinician to hand mould it to the child's foot after using a heat gun on the underside of the orthotic. It can also be heated in an oven and moulded to a cast.

**SIZING:** A - H, see page K2

Product #	Description
204	3/4 length
205	Full length
218	Narrow



204



218



## 10761 Multiflex Reinforced Insole

An orthotic blank that is heat mouldable and vacuum formable. These orthotics have a multi-form base for cushioning, a reinforcing layer for support, and a Purofit top cover for long-lasting conformity to the user's foot. Ideal for arthritic patients.

Sizes 1-7

**SIZING:** See page K2



## 10898 Plastoform Diabetic Insole

The Plastoform tri-density orthotic is a prefabricated orthotic that may be heated and vacuum formed. The top layer is plastoform—a soft foam similar to plastazote with more long-term resiliency and rebound. The middle poron layer provides additional cushioning while the bottom, reinforced multiform layer provides a solid base of support for the orthotic. Additional base material allows postings to be ground directly into the orthotic.

Heat mouldable and vacuum formable. Multiform base with dual density plastazote/poron cover. Ideal for diabetic patients.

Sizes 35-48

**SIZING:** See page K2



## HL Adjust-A-Heel Lift

This adjustable heel lift now makes correcting a patient's leg length discrepancy easier than ever! Made of three separate 1/8" layers, it can easily be adjusted in height from 3/8" to 1/4" to 1/8", simply by removing one or two of the peel off layers. This heel lift can also be used for other heel symptoms such as Achilles Tendonitis or Painful Heel Syndrome. Its compact design makes for easy placement, while its soft suede cover and firm resilient rubber base are excellent shock absorbers. Sold individually.

Size	Product #	Shoe Size	
		Men's	Women's
S	HL1		4 - 7
M	HL2	6 - 8	8 - 10
L	HL3	9 +	11 +

### 303 Intercept Heel Cushion

This gel heel cushion was designed to reduce heel pain and Achilles problems. The pad absorbs shock in the heel to reduce irritation in the heel and Achilles tendon. It also has some raise to lift the heel and shorten the Achilles tendon.

**SIZING:** Shoe size

Size	Men	Women
S	5 - 6	4.5 - 8.5
M	6.5 - 10.5	9 - 12
L	11 +	




### 304 Intercept Heel Spur Pad

Designed for varying degrees of heel spurs. Can be used for minor to moderate heel spurs. The raised edges (on bottom) support the heel and allow normal walking as it lifts the foot off of the spur. The hollowing also cushions the calcaneus at the spot of the spur.

**SIZING:** Shoe size

Size	Men	Women
S	5 - 6	4.5 - 8.5
M	6.5 - 10.5	9 - 12
L	11 +	



### 304C Intercept Heel Spur Pad with Curl

Same as the 304 Heel Spur Pad but with a "curl" to allow the weight and pressure of the foot to release, relieving direct pressure off the central heel spur or bruised area by absorbing and dispersing the shock energy.

**SIZING:** Shoe size

Size	Men	Women
S	5 - 6	4.5 - 8.5
M	6.5 - 10.5	9 - 12
L	11 +	



### 2673 DynaGel™ Comfort Heel Cup

DynaGel™ products were developed in a bioengineering laboratory to match the elasticity of human skin. This unique formulation optimizes pain relief caused by friction and pressure on feet, toes and fingers. DynaGel™ products are infused with mineral oil and skin care ingredients authenticated by the United States Pharmacopeia for skin moisturizing and healing. Antimicrobial ingredients in the gel and fabric prevent bacterial growth.

**Features:**

- Absorbs the shock normally dissipated by the natural fat pad on the bottom of the heel, making it great for running or standing on hard surfaces.
- The soft gel centre promotes the healing of a calcaneal spur.
- Comes in pairs.

**SIZING:**

Size	Length x Width
S/M	10.5 x 6.25 cm
L/XL	11.5 x 7.5 cm



## Formthotics™



Formthotics™ are a foot orthotic system made of a closed cell, non absorbing, hypo-allergenic polyethylene foam that will dry within minutes of use. These orthotics are heat mouldable, taking no more than 5-10 minutes for the complete custom process. The design features of Formthotics™ incorporate all the elements of biomechanical theory. Formthotics are moulded to the patient's feet and shoes while maintaining the functional neutral posture determined by your examination.

**\*Each style comes in select sizes, please contact us to confirm if your requested style is available in XXS/XXL/XXXL**



Teal/Charcoal



Lime/Charcoal



Purple



Red



Yellow/Charcoal



Yellow/Red



Orange/Charcoal



Green/Blue

### 306 Original Sport

The original Formthotics insole, tested and trusted across the years to provide custom support. Available in youth sizing.

Product #	Description
306	Dual (Teal/Charcoal)
306S	Single (Purple)
306Y	Youth Dual (Teal/Charcoal)

### 307 Sport Cycle

Minimalist construction specifically designed for tight fitting cycle shoes. Functional support without bulk for road, off-road, triathlon, track, and recreational cycling.

Product #	Description
307	Dual (Lime/Charcoal)
307S	Single (Purple)

### 308/9 Sport Ski

Optimum fit and support for the advanced technical or competitive skier, in tight fitting race boots. Dual is best for the more aggressive or demanding skier who wants support.

Product #	Description
308S	Single (Purple)
308H	Single High (Red)
309H	Dual High (Yellow/Charcoal)
309U	Dual Ultra (Purple)

### 308/9T Sport Ski Touring

Slim footbed adding comfort, support, and stability when walking or skiing on any terrain. The best option when the boot is a tighter fit and there is limited space available.

Product #	Description
308T	Single (Purple)
309T	Dual High (Lime/Charcoal)

### 310 Sport Football

Shaped to fit in narrower fitting sports boots. Increases comfort and protects the feet from stud pressure where space is limited. Great in skates, rugby boots, and slim-cut sneakers.

Product #	Description
310	Dual (Lime/Charcoal)
310S	Single (Purple)

### 311 Sport Run

Place feet in a natural supported position for greater connection to the ground. Increased shock absorbency and even distribution of body weight across the feet's surface reduces stress on pressure points.

Product #	Description
311	Dual (Teal/Charcoal)
311S	Single (Purple)
311C	Competition Dual (T/C)

### 312 Everyday Industry Workfit

Reliable in tough working conditions. Minimises the impact of hard surfaces and provide custom support in work-boots, spreading the load over the surface of the foot. (Orange/Charcoal)

### 313 Sport Hike Dual

Heel and arch support for best boot fit, cushion comfort, and stability on steep and craggy ground. (Teal/Charcoal)

### 316 ShockStop Diabetic Insole

Shock-absorbing designed runners, walkers, and people who spend long periods standing on hard surfaces. Beneficial for patients requiring extra care: history of arch blistering, elderly, diabetes, arthritis, and atrophy. (Green/Blue)

### 314 Original Dual Medium

The original Formthotics insole, tested and trusted across the years to provide custom support. Available in youth sizing.

Product #	Description
314	Original (Red/Blue)
314W	Wide-Fit (Green/Green)
314T	3/4 Length (Red/Blue)
314S	Slim-Fit Perforated (G/G)
314P	Perforated (Tan/Black)
314Y	Youth (Red/Blue)



Red/Blue

### 317 Original Single Layer

The original single layer Formthotics insole, made for when there is limited space in the shoe with seamless edges and a tapered toe.

Product #	Description
317W	Extra Soft (Beige)
317S	Soft (Burgundy)
317B	Medium (Blue)
317F	Firm (Black)
317R	Hard (Red)
317C	Cycle Firm (Black)
317WY	Youth Extra Soft (Beige)
317BY	Youth Medium (Blue)



Green/Green



Tan/Black

### 326 Low Volume Single

Low arch profile while still maintaining support, available in firm or with a perforated toe - along with the tapered style.

Product #	Description
326P	Medium Perforated (Green)
326F	Firm (Black)



Beige

### 327 Low Profile Dual

Ideal for patients with flat feet. Softer top layer improves comfort and reduces movement in the shoe.

Product #	Description
327	Dual (Green/Green)
327R	Dual Hard (Red/Red)



Blue

### 339 Original Dual Hard

High density semi-rigid base with a softer top layer, enhanced biomechanical support and durability with enhanced comfort and sensory input from the soft top layer.



Burgundy

### 341 Comfort Single, 3/4 Length

Gentle support helps reduce foot and leg fatigue. Top layer of foam helps prevent foot-slippage inside the boot, full contact improves balance, and the closed-cell foam means they will dry out quickly.

Product #	Description
341	Medium (Black)
341R	Hard (Red)



Black

#### DENSITY GUIDE

**Single:** The best option when space in the shoe is limited. Contours to the shape of the foot without the extra bulk.

**Dual:** An additional cushioning foam layer on top of the single's foam for added comfort. Best when space is not limited.

**High:** High density Formax foam for minimal bulk and high impact use.

**Ultra:** High density Formax foam base with additional cushioning foam layer on top. For high intensity activities.

#### SIZING: Shoe Size

Size	Men's	Women's	Size	Youth
XXS	4-5	5	Y1	13-1.5
XS	5.5-6	5.5-7	Y2	2-3.5
S	6.5-7.5	7.5-9	Y3	4-5
M	8-9	9.5-10.5		
L	9.5-10.5	11-12		
XL	11-12.5	-		
XXL	13-15	-		
XXXL	16-17	-		

## Formthotics - 3D Additions



### 328 Heel Raise

Alter biomechanics and function in the sagittal plane at the rearfoot

- May be attached under the heel area of Formthotics™ to elevate the calcaneus
- 3D cup around the heel
- Provides approximately a 8 mm raise
- Helps alleviate Achilles tendon and triceps surae strain
- Assists in balancing leg length discrepancy
- Helps improve function in a forefoot equinus
- Self-adhesive



### 321 Tri-Plane Rearfoot Wedge

Alter foot biomechanics and the triplane action of dorsiflexion, eversion and abduction at the rearfoot.

- Attaches to the under surface of Formthotics™
- May be used as a medial or a lateral wedge
- 3D cup around the heel
- Quick easy alignment with Formthotics™ “teardrop”
- Designed for accurate integration with arch support
- Tri-plane alignment of wedge
- Minimal bulk
- Self-adhesive



### 330 Longitudinal Arch Pads

Provides additional support of the medial column of the foot.

- Attaches under the arch of Formthotics™
- Easy and obvious integration with Tri-Plane Wedge and Formthotics™ “teardrop”
- Increases midfoot inversion and stabilises excessive midfoot pronation
- Utilises a softer compression-moulded EVA foam
- 3D fit to medial border
- Minimal bulk
- Self-adhesive



### 331 Reversible Forefoot Wedges

An orthotic addition suitable for forefoot modification of Formthotics™.

- Can be reversed to use on either lateral or medial aspect
- Easy to modify
- Long, gentle taper
- Self-adhesive

## Formthotics - Orthotic Components



### 318 Heel Raise

Available in 4mm or 6mm thickness.  
1 or 5 pairs/pack  
Sizes: S-L




### 319 Arch Pad

1 or 5 pairs/pack  
Sizes: S-L



### 322 Heel Wedge

Available in firm and medium density.  
1 or 5 pairs/pack  
Sizes: S-L



### 323 Extended Wedge

Available in firm and medium density.  
1 or 5 pairs/pack  
Sizes: S-L



### 324 Metatarsal Dome

Available in firm and medium density.  
5 pairs/pack  
Sizes: 1 and 2

## Formthotics - Orthotic Components

### 338 Wedged Strip

Designed to streamline Formthotics customization for clinical efficiency and patient satisfaction.

Offers easy, quick, full-length modifications, with minimal mess and waste.

- Made from Formthotics medium PE140 Formax foam with self-adhesive backing
- 300mm length, 45mm width, 4mm height with a 5 degree angle
- Attaches easily to the lateral or medial underside of Formthotics
- Helps increase support at the edges to counteract forces affecting the patient's gait
- Available in packs of 3 pairs, a total of 6 units per pack



---

### FHMV20200 Formthotics Heating Machine

The best tool to use for custom moulding. It is designed to provide the perfect fit in just 3 minutes. Supports achieving the optimal result with Formthotics thermo-formable polyethylene foam to be moulded to the shape of the individual's foot and shoe.



### FHMV20BE Formthotics Heater Boot Extension

Designed for use with ski boots, snowboard boots, or other fully sealed footwear.

Used with the Formthotics Heating Machine, it gently heats any orthoses to the optimal temperature for custom moulding to the foot and footwear.



## Tuli's Footcare Products

**TULI'S®**

Designed to duplicate nature's own shock absorbing system through the use of a waffle heel cup and shock-absorbing material, eliminating the shock of walking and running with each step. Tuli's heel cups can relieve heel spurs, shin splints, sore aching feet, painful legs and knees, arthritic joints, and back and neck pains.



### 346/347 Tuli's Gel Heel Cup

A gel polymer designed to provide a softer waffle heel cup. Can be worn in dress or casual shoes for all activities. Provides fast relief for a painful heel problem. Retains absorbency/cushioning effect. Will not bottom out. Offers protection from heel pain and shock. Use in athletic or work shoes. Fully supports the body in all sports and work related activities. The Heavy Duty Heel Cup (347) offers higher and heavier ribs than standard models. Helps improve performance at work and play.

Product #	Description
346	Regular Gel Heel Cup
347	Heavy Duty Gel Heel Cup

**SIZING:**

Size	Weight
S	< 125 lbs (#347 only)
Reg	125-175 lbs
L	over 175 lbs



### 351 Tuli's Classic Heel Cup

Worn in your shoes, Tuli's will help prevent sore heels and knees, shin splints and heel spurs. They will also relieve arch, ankle, leg, and back pain that result from everyday heel strike. Double-ribbed design for extra shock absorbing protection. Tuli's Pro Heel Cups offer instant relief from: tired, burning aching feet, and legs, shin splints, painful heels, heel spurs, arthritic joints, painful ankle joints, nagging backache, and swollen ankles.

**SIZING:**

Size	Weight
Reg	175 lbs and under
L	over 175 lbs



### 358 Tuli's Heavy Duty Heel Cup

Tuli's Heavy Duty Heel Cups™ are twice as thick as Tuli's® Classic Heel Cups™ for maximum support and protection. They deliver the heavy duty protection you need from heel strike pain during sports activities such as volleyball, soccer, basketball and tennis. These heel cups work best in athletic shoes, cleats, work boots and more spacious lace up shoes.

**SIZING:**

Size	Weight
S	less than 80 lbs
Reg	80-175 lbs
L	over 175 lbs



### 352 Tuli's So Soft Heel Cup

Formulated to be lighter, softer, and more resilient than other heel cups. These new heel cups have a soft cloth fabric cover for ease of sliding the foot into the shoe and more comfort with or without socks. The Tuli's multi-cell design absorbs shock just like the system naturally found in your feet. This shock absorption and elevation of the heel bone provides immediate pain relief and assists in long-term healing. The low-profile design allows them to be worn in almost any shoe.

**SIZING:**

Size	Weight
Reg	up to 175 lbs
L	over 175 lbs

### 359 Tuli's Cheetah

Lightweight, breathable, Neoprene ankle support for comfort and stability. Built-in shock absorbing Tuli's Heel Cup to absorb impact energy. Easy to wear, gentle and comfortable on your feet. These heel cups are a must for gymnasts or any sports-minded person. Washable by hand or machine. Comes in beige. Universal size comes with straps, XS - XL is a sleeve.

**SIZING:** Instep circumference (around the arch of the foot)

Size	Instep Circumference	Size	Instep Circumference
XS	up to 7.5 inches	L	9.5" - 10.5"
S	7.5" - 8.5"	XL	10.5"+
M	8.5" - 9.5"	U	Universal



### 357 Tuli's Metatarsal Cushion

Provides instant relief for tired, swollen, aching metatarsals. Tuli's metatarsal cushions relieve that agonizing forefoot pain when wearing high heel fashion or dress shoes.

**SIZING:** Universal



### 349 Tuli's Plantar Fasciitis Insole

The Tuli's Plantar Fasciitis Insoles are full-length premium arch supports designed specifically for plantar fasciitis. The arch supports and the metatarsal arches evenly distribute body weight pressure throughout the foot to help improve stability, balance and posture. A 4° angled medial heel wedge is incorporated to encourage stability by preventing excessive over-pronation, which can cause stress and inflammation on the plantar fascia. The cushion in the heel and forefoot incorporates the patented Tuli's multi-cell shock absorption technology for natural shock absorption, providing long term wear comfort. Fits in any work boot, running, dress and casual shoe. The front of the insole can be trimmed for optimum fit.

**SIZING:** Shoe size

Size	Men	Women
S	5.5 - 7	6.5 - 8
M	7.5 - 9	8.5 - 10
L	9.5 - 11	10.5 - 12
XL	11.5 - 13	12+





**A**



**B**



**C**




**D**



**E**



**F**



**G**

**H**



**I**

**(A) 5580 Hallux Valgus Night Splint**

The Hallux Valgus night splint is designed to gently push your toe back to the natural position to help relieve pain over time. This is a rigid splint which should not be worn during weight-bearing activities.

**SIZING:** Universal, specify left or right

**(B) 479 Leather Metatarsal Pad - Kidney**

Soft foam rubber metatarsal pad with leather top cover and an adhesive backing for easy application. Kidney shaped. Pairs.

**SIZING:** 1 - 4

**(C) 20770 Leather Metatarsal Pad - Teardrop**

Soft foam rubber metatarsal pad with leather top cover and an adhesive backing for easy application. Tear drop shape. Pairs.

**SIZING:** 1 - 4

**(D) 112070 Metatarsal Bandage**

**SIZING:** Forefoot circumference 16 - 26 cm

**(E) 810 Hallux Valgus Daytime Splint**

Elastic Hallux Valgus support that will fit inside the shoe for daytime wear. Specify left or right.

**SIZING:** S, M, L

**(F) 26140 Hallux Valgus Night Splint**

High quality splint for nighttime correction of Hallux Valgus.

Size	Product #	
	Left	Right
Child	<a href="#">26140-341</a>	<a href="#">26140-340</a>
Women	<a href="#">26140-371</a>	<a href="#">26140-370</a>
Men	<a href="#">26140-411</a>	<a href="#">26140-410</a>

**(G) 823 Toe Crest Cushion**

A comfortable fitting foam toe spacer.

**SIZING:** S, L

**(H) 829 Buddy Toe Splint Cushion**

Foam, plush, Velcro® sensitive toe splint. To splint one toe to the next or provide padding between the toes.

**SIZING:** Universal

**(I) 180 Hallux Valgus Splint**

The **180** is a soft, comfortable daytime Hallux Valgus support that can be worn in the shoe. Made of high-grip elastic Neoprene, it will provide a gentle stretch all day and night for the big toe in the valgus position. Specify left or right.

**SIZING:** Universal



## Luga COMFORSIL Footcare Products

Luga items have been discontinued, limited stock available. Please contact to confirm.

All Comforsil products contain an antiseptic solution produced from the product Calbenium® which is anti-bacterial, antiviral and antifungal. All Luga products come in retail packages.

### CC206 Toe Separators



Provide comfort and relief for painful and annoying soft corns found between the toes. Relieves discomfort caused by overlapping toes. The agreeable density and properties of high viscosity Comforsil provide comfort and protect the sensitive toes. 2 per pack.

**SIZING:** S, M, L

### CC218-2 Toe Spreader



A corrective shield between first and second toes, the extra soft skin coloured Comforsil toe spreader effectively prevents deviation of the hallux into a resulting bunion, while avoiding friction. Ideal in rehabilitating the first toe after a unionectomy. 2 per pack.

**SIZING:** S, M, L

### CC229 Elastic Band with Silicone Metatarsal Insert



This product follows the design of the more traditional elastic metatarsal band to address met head problems, both biomechanical, and/or lesions. An extra soft Comforsil disc has been embedded in the plantar aspect of the elastic band to help painful metatarsals unload and diffuse pressure. It fits snugly in any kind of shoe. Single.

**SIZING:** S to XL, forefoot circumference

### CC256 Forefoot Cushions



Ideal for the protection and pain relief of forefoot conditions. They help reduce calluses as they afford a therapeutic and relaxing effect on the foot. Made of soft elastic fibre. Place tab between first and second toes. Washable. Simply air dry the product after washing. One pair.

**SIZING:** Small, Large

### CC220 Metatarsal Pad with Reinforced Loop



This extra soft, skin coloured Comforsil metatarsal pad is intended to help eradicate plantar calluses and hard skin, reduces the adverse effects of weight bearing, and alleviates painful rubbing and friction when walking. Specifically designed to treat forefoot soreness, and as a preventive device for those who stand long hours on their feet. An elastic fabric, embedded in the digital loop, reinforces it for added durability. One pair.

**SIZING:** Universal

### CC226 Toe Crest with Reinforced Loop



A variation on the silicone toe prop, the toe crest version presents with a reinforced digital loop, and addresses hammer toe conditions. The reinforcement of the loop consists of an elastic fabric woven within the silicone for added strength and durability. As with all splints that incorporate a loop, we emphasize the need for handling delicately. Single.

**SIZING:** S, M, L

### CC215 Heel Spur Cushions with Removable Cavity



The convenient, removable dot of this Comforsil Heel Spur Cushion permits the placing of the painful spur directly into the Comforsil cut out, thus completely unloading the troublesome calcaneal spur, especially in very acute cases. The device also distributes body weight away from the affected area. One pair.

**SIZING:** S, M, L

### CC227-2 Toe Caps



A variation on the silicone toe prop, the toe crest version presents with a reinforced digital loop, and addresses hammer toe conditions. The reinforcement of the loop consists of an elastic fabric woven within the silicone for added strength and durability. As with all splints that incorporate a loop, we emphasize the need for handling delicately. Single.

**SIZING:** S, L

### CC212 Heel Spur Cup



A preventive, dual density, extra soft, skin coloured Comforsil heel cup for the tender hind foot. The base of the calcaneus is cradled within an extra soft circular spot of Comforsil silicone. One pair.

**SIZING:** S



20781



20782

## Novaped Adjustment Wedges

For individual adjustment of Novaped foot supports, these elements can be quickly, securely and durably applied without any adhesive product. The Pronation/Supination Wedge can be applied to the medial or lateral side of the heel or forefoot for a quick and easy posting option. The Heel Wedge provides a quick solution for heel lift or posting.

Size 1 for insole size 35-39 and size 2 for insole size 40-43.

Product #	Description
20781	Pronation/Supination Wedge 4mm
20782	Heel Wedge 4mm



22400



22500



22600

## PPT Self Adhesive Components

Self Adhesive PPT® Components allow for quick and easy application. (1 set includes 6 pairs)

**Longitudinal Arches** - Provides longitudinal arch support. Recommended for flat insoles with flange.

Product #	Description
22400	PPT Longitudinal Arch Pads 6 prs/pk (Small)
22410	PPT Longitudinal Arch Pads 6 prs/pk (Medium)
22420	PPT Longitudinal Arch Pads 6 prs/pk (Large)

**Metatarsal Bars** - Relieves pressure across the metatarsal heads.

Product #	Description
22500	PPT Metatarsal Bars 6 prs/pk (Small)
22510	PPT Metatarsal Bars 6 prs/pk (Medium)
22520	PPT Metatarsal Bars 6 prs/pk (Large)

**Metatarsal Pads** - Transfers forces from the metatarsal heads.

Product #	Description
22600	PPT Metatarsal Pads 6prs/pk (Small)
22610	PPT Metatarsal Pads 6prs/pk (Medium)
22620	PPT Metatarsal Pads 6prs/pk (Large)



### 31150-31155 Toe/Finger Spreader

**This item has been discontinued, limited stock available. Please contact to confirm.**

A soft silicone pad that is molded to fit comfortably between toes or fingers. Anti-bacterial, hypoallergenic, and non-toxic. Facilitates proper toe or finger alignment and prevents overlapping digits.

Product #	Size	Product #	Size
31150	Small - 4/bag	31153	Small - 50/bag
31151	Medium - 4/bag	31154	Medium - 50/bag
31152	Large - 4/bag	31155	Large - 50/bag



### 31160-31163 Toe Separator

**This item has been discontinued, limited stock available. Please contact to confirm.**

Protects against interdigital rubbing and irritation. Anti-bacterial, hypoallergenic, and non-toxic. Comforts between toe or finger formations such as corns or warts.

Maintains alignment during healing for extended and post-operative use. Place separator so the narrowest part is between the toes/fingers and towards the foot/hand.

Product #	Size
31160	Small - 15/bag
31161	Large - 15/bag

## Intelligel Performance Achilles Protector



**This item has been discontinued, limited stock available. Please contact to confirm.**

This lightweight, comfortable, toeless sleeve does not create bulk in the shoe and can be worn with casual, dress, and work shoes. Its anatomically-shaped soft gel pad comforts, cushions, and protects the tender Achilles area. Recommended by medical professionals for people with Achilles tendonitis, bursitis, heel prominences, pump bumps and hardened or cracked skin. Washable and reusable.

Product #	Size
35038	S/M
35039	L/XL



## Intelligel Performance Boot Bumper



**This item has been discontinued, limited stock available. Please contact to confirm.**

Conforming 1/8" polymer gel pads protect the front and back of the ankle from pressure, pinching and friction. Irritation, burning and redness are soothed by the unique properties of the Intelligel. Elastic sleeve protects the foot and ankle from abrasions caused by athletic footwear without the added bulk. Ideal for use in ice skates or ski boots. Comes in retail packaging.

Product #	Size
35186	S
35187	L



## 35084 Intelligel Performance Shin Guard



**This item has been discontinued, limited stock available. Please contact to confirm.**

This 3" x 5" elastic tube is lined with gel to reduce pressure, pinching, or friction on the lower leg. Ideal for use in ice skates or ski boots. One per package. Comes in retail packaging.

**SIZING:** Universal



## M-Gel Body Disc



**This item has been discontinued, limited stock available. Please contact to confirm.**

This unique body disc allows for a totally new approach for the prevention/treatment of pressure and friction induced bruises and lesions. Medical grade mineral oil provides comfort. 2 discs per package.

Product #	Description
31030	2.5" diameter
31031	4" diameter



## DynaGel™

Ortho Active is proud to introduce our new footcare line: DynaGel™. DynaGel™ products were developed in a bioengineering laboratory to match the elasticity of human skin. This unique formulation optimizes pain relief caused by friction and pressure on feet, toes and fingers. DynaGel™ products are infused with mineral oil and skin care ingredients authenticated by the United States Pharmacopeia for skin moisturizing and healing. Antimicrobial ingredients in the gel and fabric prevent bacterial growth.

The mineral-oil-infused gel protects and comforts the toes. Medical grade mineral oil is slowly released to soften and moisturize the skin. Long lasting gel will not flatten or wear out like foam products.



### 2600 Digital Tube

- Mineral-oil-infused antibacterial gel provides padding, relieves pressure and eliminates friction. Medical grade mineral oil is slowly released to soften the skin, helping to relieve scar tissue.
- Hyper-elastic fabric is conformable and comfortable for any shape.

**2600-5** - Five 15cm tubes per bag

**R2600-4** - Four 2cm tubes per retail box

Size	Diameter
S	0.5 cm
M	1 cm
L	1.5 cm



### 2626 Digit Tube

- Antibacterial gel provides padding, relieves pressure and eliminates friction.
- This soft, extra-durable gel does not need a fabric backing; it stretches easily and provides less compression than our regular ribbed tube.
- Easily cut to fit required length.

**2626-10** - Ten 3cm tubes per bag

**R2626-2** - Two 3cm tubes per retail box

Size	Diameter
S	0.5 cm
M	0.8 cm
L	1 cm



### 2647 Hammer Toe Support

- The DynaGel Hammer Toe Support is designed to help straighten bent toes and cushion toe top corns caused by bent toes.
- DynaGel will soften and smooth corns, relieve pain from bent-under toes, and make walking easier and more comfortable.

**2647-10** - Ten per bag

**R2647-2** - Two per retail box

Size	Diameter
S	0.6 cm
M	0.8 cm
L	1 cm



### 2601 Digital Cap

- Mineral-oil-infused antibacterial gel provides padding, relieves pressure and eliminates friction. Medical grade mineral oil is slowly released to soften the skin, helping to relieve scar tissue.
- Hyper-elastic fabric is conformable and comfortable for any shape.
- Easily cut to fit required length.

**2601-10** - Ten 6cm caps per bag

**R2601-2** - Two 6cm caps per retail box

Size	Diameter
S	0.6 cm
M	0.8 cm
L	1.5 cm



### 2633 Digital Cap

**This item has been discontinued, limited stock available. Please contact to confirm.**

- Antibacterial gel provides padding, relieves pressure and eliminates friction. This soft, extra-durable gel digital cap does not need a fabric backing, it stretches easily, and provides less compression than our regular ribbed cap.
- Easily cut to fit required length.

Size	Diameter
S	0.5 cm
M	1 cm
L	1.5 cm

**2633-10** - Ten 4cm caps per bag

**R2633-2** - Two 4cm caps per retail box



### 2637 Toe Spreader with Loop

- Gel toe spreader re-aligns toes and relieves pressure from bunions.
- The loop keeps the toe spreader in the correct position, helping to prevent the toes from overlapping.

**2637-10** - Ten per bag

**R2637-2** - Two per retail box

Size: Universal



## 2608 Digital Pad

- An easy-to-cut elastic tube for making individual pads when buying on a strip
- Mineral-oil-infused antibacterial gel provides padding, relieves pressure and eliminates friction.
- Hyper-elastic fabric is conformable and comfortable for any shape, and holds the pad in place to wear easily on fingers and toes.

Per bag:

Sizes: S, M, L

**2608S-16** - Two strips, 8 pads per strip (pad diameter: 1cm)

**2608M-14** - Two strips, 7 1.5cm pads per strip (pad diameter: 1.5cm)

**2608L-12** - Two strips, 6 2cm pads per strip (pad diameter: 2cm)

**R2608-2** - Available in S/M/L, two individual pads per retail box



## 2628 Toe Spreader with Two Loops

**This item has been discontinued, limited stock available. Please contact to confirm.**

- Gel toe spreader re-aligns toes and relieves pressure from bunions.
- The two loops keep the toe spreader in the correct position, helping to prevent the toes from overlapping.

Size: Universal

**2628-10** - Ten per bag

**R2628-2** - Two, retail box



## 2631 Five Toe Separator

- Gel toe separator re-aligns toes and relieves pressure between overlapping toes where more than two toes on the same foot are involved.

**2631-10** - Ten per bag

**R2631-2** - Two per retail box

Size: Universal



## 2604 Toe Separator

- Gel toe separator re-aligns toes and relieves pressure for under/overlapping toes.
- The anatomical shape allows the separator to fit easily in footwear to minimize toe pressure.

**2604-15** - 15 per bag / **2604-50** - 50 per bag

**R2604-2** - Two per retail box

Size: S, M, L



## 2605 Toe Spreader

- Gel toe spreader re-aligns toes and relieves pressure from bunions.

**2605-15** - Fifteen per bag

Sizes: S, M, L

**2605-50** - Fifty per bag

**R2605-2** - Two per retail box



## 2607 Bunion Sleeve

**This item has been discontinued, limited stock available. Please contact to confirm.**

- The conformable gel pad of this Bunion Sleeve absorbs shear forces and cushions the bunion, while helping to realign the big toe.
- The thin, sheer fabric will hold the bunion pad in the correct position and fits easily in any shoe. Can be worn on either left or right foot.

**R2607** - One per retail box

Size	Width
S/M	8 cm
L/XL	9 cm



## 2615 Toe Spreader w/ Sleeve

**This item has been discontinued, limited stock available. Please contact to confirm.**

- Gel toe spreader re-aligns toes for Hallux Valgus correction and relieves pressure.
- The ultra-thin and sheer elastic sleeve holds the spreader in the correct position, and will fit easily in any shoe.

One per retail box

Sizes: S/M, L/XL

See **2607** for sizing



## 2609 Corn Pad

- Hyper-elastic fabric is conformable and comfortable for any shape.
- Strips can be easily cut to fit any length requirement.

**2609-3** - Three 25cm strips Sizes: S, M, L

**R2609-3** - Three 3cm pads See **2647** for sizing



## 2630 Bunion Guard

**This item has been discontinued, limited stock available. Please contact to confirm.**

- This anatomically-shaped, mineral-oil-infused gel bunion pad relieves pressure, absorbs shear forces and cushions the bunion for comfort and pain relief. It also helps realign overlapping toes.
- The toe loop helps to keep the bunion guard in the correct position at all times. Can be worn on left or right foot.

**R2630-2** - Two per retail box Size: Universal



## 2629 Bunion Guard with Toe Spreader

- Gel toe spreader re-aligns toes to prevent them from overlapping for Hallux Valgus correction, relieving pressure from bunions. By attaching the spreader to the bunion guard, the spreader is always held in the correct position. Can be worn on left or right foot.
- The anatomically shaped gel bunion pad absorbs shear forces and cushions the bunion for comfort and pain relief.

**R2629-2** - Two per retail box Size: Universal



## 2634 5<sup>th</sup> Digit Valgus Guard

**This item has been discontinued, limited stock available. Please contact to confirm.**

- This gel toe spreader realigns toes to prevent them from overlapping, providing 5th digit Hallux Valgus correction and relieving pressure from 5th digit bunions. The anatomically shaped gel bunion pad will absorb shear forces and cushion the bunion. Can be worn on left or right foot.
- The spreader is attached to the bunion guard, so the spreader is always held in the correct position.

**2634-10** - Ten per bag Size: Universal  
**R2634-2** - Two per retail box



## 2617 Toe Spreader/Bunion Protector

**This item has been discontinued, limited stock available. Please contact to confirm.**

- The gel pad on this Bunion Sleeve absorbs shear forces and cushions the bunion.
- Gel toe spreader re-aligns toes for Hallux Valgus correction and relieves pressure
- The ultra-thin, sheer elastic sleeve holds the spreader in the correct position and fits easily in any shoe. Can be worn on left or right foot.

**2617-10** - Ten per bag Sizes: S/M, L/XL  
**R2617** - One per retail box See **2607** for sizing



## 2632 Toe Crest

- This soft gel toe crest provides effective cushioning to relieve the discomfort of crooked toes, hammer toes, mallet toes and claw toes.
- The elastic toe loop is designed for wearing comfortably.
- Slide the loop over the second or third toe for best effect and comfort.

**2632-5** - Five pairs per bag Size: Universal  
**R2632** - One pair per retail box



## 2636 Metatarsal Pad with Loop

- The soft gel of this DynaGel metatarsal pad absorbs shock and reduces pressure to relieve metatarsal pain and discomfort.
- The elastic, thin toe loop keeps the pad in place and is easy to apply.
- Long lasting gel will not flatten or wear out like foam products.

**R2636-2** - Two per retail box Size: Universal



## 2622 Heel/Elbow Protector

**This item has been discontinued, limited stock available. Please contact to confirm.**

- The comfortable gel pad incorporated into the back of our DynaGel sleeve reduces pressure and friction where the Achilles attaches to the heel to alleviate this painful condition.
- The anatomically-shaped open-toe elastic sleeve is designed to hold the gel pad in the correct position.
- This sleeve will fit most ankle sizes, left/right.

One per retail box Size: Universal

### 801 Sportec Compression Socks

The Sportec Compression Socks have a maximum compression level of 30mmHg, and have been designed with several features for high performance activity. Targeted achilles and plantar fasciitis compressions zones provide extra support of the achilles and foot arch. Heel and toe Nano-Silver antimicrobial fibre provide odour reduction, and the loose-knit toe prevents blisters and abrasion. Made with moisture-wicking construction for cool performance.

**SIZING:** (in cm)



Size	S	M	L	XL
cC	28-38	33-43.5	35.5-46	38.5-51.5
cB	19-21.5	21.5-24	23-26.5	26.5-32



### 5192 FXT Compression Socks



Designed with targeted compression zones to aid in the pain management and relief from plantar fasciitis, heel spurs and arch pain. Ideal for sports or those on their feet all day.

**SIZING:**

Shoe Size	XS	S	M	L	XL
Men's	4 - 5.5	6 - 7.5	8 - 9.5	10 - 11.5	12 - 14
Women's	5 - 6.5	7 - 8.5	9 - 10.5	11 - 12.5	13 - 15



### 70U Soft Plantar Fasciitis Night Splint

For nighttime treatment of plantar fasciitis, this Ortho Active product is designed to be worn while sleeping to provide a soft stretching action and alleviate foot pain caused by plantar fasciitis upon waking. Included elastic straps provide a long-term, low-load stretch, and the anti-slip sole makes it easy for patients to wear the splint around the house. Fits left or right.

**SIZING:** Universal





### 6070 PF Night Sock

This is a comfortable, low-profile support made with lightweight and breathable materials. The Night Sock helps relieve the heel pain associated with common overuse injuries such as Plantar Fasciitis and Achilles Tendonitis when taking those first morning steps. Has a comfortable knitted sleeve design and an adjustable dorsiflexion strap which gently pulls the toes to stretch the calf muscle and soft tissues of the foot and ankle.

**SIZING:**

Size	Calf Circumference
Regular	Up to 17"
Large	18" - 22"



### 73 Sole Stretch

The Ortho Active Sole Stretch provides targeted stretching to alleviate plantar fasciitis during a low-load prolonged stretch. The amount of stretch for the ankle and metatarsal joints is independently adjustable to suit each person's needs. The soft design provides additional comfort while at rest.

**Indications:**

- Plantar Fasciitis
- Plantar flexor tightness or contracture

**SIZING:** Shoe size, calf circumference 1" below the knee

Size	Women's Shoe Size	Men's Shoe Size	Calf Circ.
S	5-8	6	9"-15"
S-CM	5-8	6	11"-17"
M	9-11	7-10	11"-17"
M-CL	9-11	7-10	13"-19"
L	12-14	11-14	13"-19"
L-CXL	12-14	11-14	17"-24"



### 5380 Plantar Fasciitis Adjustable Night Splint

Plastic night splint with adjustable dorsiflexion angle, posterior cutouts for breathability and comfort, and non-skid rubber sole material. Removable, washable, breathable, 1/2" thick foam liner with soft Velcro®-sensitive liner cover. Removable toe elevation pad. Padded shin and ankle straps with quick release buckles. Provides a low-load prolonged stretch to the plantar fascia and achilles tendon.

**Indications:**

- Plantar Fasciitis
- Achilles Tendonitis
- Plantar flexion contractures

**SIZING:** Shoe size

Size	Men	Women
S	5 - 7	5.5 - 7.5
M	7 - 9	7.5 - 9.5
L	9 - 12	9.5 - 12.5



### 6576U Dorsal Night Splint

The Dorsal Night Splint is designed to hold the foot in a neutral position during sleep, thereby helping to alleviate the morning pain and discomfort associated with plantar fasciitis. This comfortable brace gives increased patient compliance due to its low profile, reduced bulk and weight, and low heat retention.

**SIZING:** Universal

### 595 Dynamic Post-Op Shoe

The Dynamic Post Operation Shoe is made from a breathable material with plush lining for comfort against the patient's skin. The included toe protector is removable for access, cleaning, or comfort, and Velcro straps built in to the product provide a wide range of adjustability to fit on any foot.

**SIZING:** Shoe size

Size	Sole Length	Men's Shoe Size	Women's Shoe Size
XS	20.5-22 cm	<4.5	<5.5
S	22.5-24.5 cm	4.5-6.5	5.5-7.5
M	25-26.5 cm	7-8.5	8-9.5
L	27-28.5 cm	9-10.5	10-11.5
XL	29-30.5 cm	11-13	12-13.5



### 596 Dynamic Heel Offloading Shoe

The Dynamic Heel Offloading Shoe has all the features of the Post Operation Shoe with a heel offloading design feature to reduce pressure on the heel. This shape helps speed the healing process after surgery and reduce the stress of walking.

**SIZING:** Shoe size

Size	Sole Length	Men's Shoe Size	Women's Shoe Size
XS	20.5-22 cm	<4.5	<5.5
S	22.5-24.5 cm	4.5-6.5	5.5-7.5
M	25-26.5 cm	7-8.5	8-9.5
L	27-28.5 cm	9-10.5	10-11.5
XL	29-30.5 cm	11-13	12-13.5



### 597 Dynamic Forefoot Offloading Shoe

The Dynamic Forefoot Offloading Shoe has all the features of the Post Operation Shoe with a forefoot offloading feature to reduce forefoot pressure. This shape assists in speeding up healing after surgery and reducing the stress of walking.

**SIZING:** Shoe size

Size	Sole Length	Men's Shoe Size	Women's Shoe Size
XS	20.5-22 cm	<4.5	<5.5
S	22.5-24.5 cm	4.5-6.5	5.5-7.5
M	25-26.5 cm	7-8.5	8-9.5
L	27-28.5 cm	9-10.5	10-11.5
XL	29-30.5 cm	11-13	12-13.5





Product No.	Product	Page
-------------	---------	------

### Pediatric Supports

983192-983124	Aspen Pediatric Collar Set	J2
P1214	Arm Sling	J2
60145	Arm Sling with Padded Strap	J2
60125	Shoulder Immobilizer	J2
3957	Wrist Splint	J3
P5155	Walking Boot	J3
8050	Easy Walk KITE AFO	J3
FH201-FH202	Trulife Child's Play Matrix AFO	J4
13614-13616	Orthomerica UFO	J4

### Children's Foot Orthotics (OTS)

230	Fashmo Flat Blanks	J4
219/220	Novatherm Multi Heel Cup	J4
204, 205, 218	Supinara Children's Foot Orthotics	J5

### Ortho Active Pediatric Supports

P32	Pediatric Hinged Knee Brace	J6
P34	Pediatric Jumper's Knee Brace	J6
P97	Pediatric Contoured Wrist Brace	J6

## Pediatric Foam Collar



This cervical collar is made of breathable, high quality foam padding to immobilize and support the neck as it heals. Soft, lightweight, and easy to wear. Has a slight contour to comfortably fit under the chin. Comes in 6cm and 8cm widths.

**SIZING:**

Product # (6cm)	Neck Circumference	Product # (8cm)	Neck Circumference
P4506S	25 - 30cm	P4508S	35 - 40cm
P4506M	30 - 35cm	P4508M	40 - 45cm
P4506L	35 - 40cm	P4508L	45 - 50cm

## 1214 Pediatric Arm Sling



**This item has been discontinued, limited stock available. Please contact to confirm.**

**Features:**

- Soft poly-cotton
- Fun dinosaur print
- Adjustable shoulder strap with a side buckle, hook and loop closure allows for customized fit

**SIZING:**

Size	Forearm Length
LG	7" - 11"

## 6014S Arm Sling with Padded Strap



This durable polyester sling is ideal for immobilization and support of the shoulder and elbow joints. Fits right or left arm and helps prevent shoulder rotation. Shoulder straps are a plush foam laminate with contact closure which provides additional comfort. The canvas pocket is 7" deep for comfortable support, and an O-Ring keeps the strap centered.

**SIZING:** Fits arm length of <13" from elbow to fingertip

## 6012S Shoulder Immobilizer



Ideal for immobilization and support of the shoulder and elbow joints. This high quality, durable 100% cotton sling helps prevent shoulder rotation. Has contact closure adjustment on shoulder and body straps. Waist and shoulder straps are a plush foam, laminated with contact closure which provides additional comfort. Waist strap is easily adjusted or can be removed. The canvas pocket is 7" deep for comfortable support, and an O-Ring keeps the strap centered.

**SIZING:** Fits arm length of <13" from elbow to fingertip

## 3957 Wrist Splint

**This item has been discontinued, limited stock available. Please contact to confirm.**

Poly cotton bound to foam and lined with nylon tricot. Pre-shaped removable stays with hook/loop closure. Specify right or left.

**SIZING:**

Size	Guidelines
Infant	6 - 18 months
Pediatric	18 months to 2" palm width



## P5155 Pediatric Walking Boot

It is designed for the treatment of stable fractures and ankle sprains for kids. It is used for early rehabilitation of a variety of post surgical and post trauma indications requiring fixed control.

**Features:**

- Rocker bottom and low profile heel facilitates a more natural gait pattern
- Lightweight design includes anatomically shaped plastic uprights that help to facilitate a better fit
- Reinforced uprights and rigid side panels protect and stabilize the ankle and lower leg
- Easy to apply and adjust with a simple strapping and liner system

**SIZING:**

Size	Foot Length	Shoe Size (US)	Shoe Size (Euro)
S	Up to 6"	Up to 8.5	Up to 25
M	6" - 7"	8.5 - 11.5	25 - 29
L	7" - 8"	11.5 - 1.5 (Youth)	29 - 33



## 8050 Easy Walk KITE AFO

intended for children with dorsiflexion weakness and is specifically designed to control medial-lateral movement and stabilize the ankle. At the same time, the carbon footplate assists the gait by flexing and providing the energy return for toe-off.

**SIZING:**

Size	Shoe Size	Height
S	20 - 24 cm	21 cm
M	24 - 27 cm	25 cm
L	27 - 30 cm	29 cm
XL	30 - 33 cm	34 cm



## Trulife Child's Play Matrix AFO



Lightweight, carbon composite frame with enhanced stability. Lateral strut is anterior to malleolus, and has an open calcaneal design. The anterior shell is height adjustable. Available with camo blue or butterfly pink shell.

Left #	Right #	Size	Colour	Height	Foot Length
<a href="#">FH201L-03</a>	<a href="#">FH201R-03</a>	M	Camo Blue	26 cm / 10 in	18 cm / 7 in
<a href="#">FH201L-04</a>	<a href="#">FH201R-04</a>	L	Butterfly Pink	29 cm / 11.5 in	20 cm / 8 in
<a href="#">FH202L-03</a>	<a href="#">FH202R-03</a>	M	Camo Blue	26 cm / 10 in	18 cm / 7 in
<a href="#">FH202L-04</a>	<a href="#">FH202R-04</a>	L	Butterfly Pink	29 cm / 11.5 in	20 cm / 8 in

## Orthomerica UFO Pediatric



The pediatric UFO™ from Orthomerica is a comfortable, prefabricated orthosis for stretching and maintaining dorsiflexion range of motion. This unique lower-limb orthosis positions the foot and ankle in optimal alignment for placing stretch on the soleus while the patient is sleeping. When used in conjunction with a knee immobilizer, this orthosis can also stretch the two-joint gastrocnemius group.

**SIZING:** Shoe size

Size	Product #	Shoe Size	Height
S	<a href="#">13614</a>	0 - 2	6.5"
M	<a href="#">13615</a>	2 - 4	8"
L	<a href="#">13616</a>	4 - 7	10.5"

## 230 Fashmo Flat Blanks (Adult and Children)



These fashmo orthotic blanks have a unique shape that provides a nice heel cup and flanges when moulded. Perfect for making UCBL-style foot orthotics. The top and bottom materials create a superb finish on a custom orthotic in less time. Available in pediatric and adult sizes.

**SIZING:** Adult size 35-46 (See size chart on page J2)  
Children size A-H - shoe size, see [83](#) chart above

**Heating Info:** Heat at 375°F for approximately 3 mins. prior to vacuum forming.

## 219/220 Novatherm Multi Heel Cup



A full-length prefab orthotic that can be remoulded. High heel counters are ideal for providing calcaneal stability. Consists of a pink or black multiform top, thermaforte bottom, and reinforcing middle layer. Available in pink or black.

**SIZING:** Children's shoe size, see [83](#) chart above

Product #	Description
<a href="#">219</a>	Multi Heel Cup Pink
<a href="#">220</a>	Multi Heel Cup Black

## 204 Suprinara Children's Foot Orthotics - 3/4 Length

Suprinara construction with high lateral support and standard heel cup. 3/4 length orthotic with multiform topcover. Heat will not transfer through multiform topcover, allowing the clinician to hand mould it to the child's foot. Works great off-the-shelf, but customization may be achieved by using a heat gun on the underside of the orthotic. Item 222 comes without top cover.

**SIZING:** Children's shoe size

Size	Children's Shoe Size	Size	Children's Shoe Size
A	4 - 5	E	12 - 13
B	6 - 7	F	1 - 2
C	8 - 9	G	3 - 4
D	10 - 11	H	5 - 6



## 205 Suprinara Children's Foot Orthotics - Full Length

Suprinara construction with standard lateral support and deep heel cup. Full length orthotic with multiform topcover. Heat will not transfer through multiform topcover, allowing the clinician to hand mould it to the child's foot. Works great off-the-shelf, but customization may be achieved by using a heat gun on the underside of the orthotic.

**SIZING:** Children's shoe size

Size	Children's Shoe Size	Size	Children's Shoe Size
A	4 - 5	E	12 - 13
B	6 - 7	F	1 - 2
C	8 - 9	G	3 - 4
D	10 - 11	H	5 - 6



## 218 Suprinara Children's Foot Orthotics - Narrow

Suprinara construction with standard lateral support and deep heel cup. Full length, narrow fit orthotic with multiform top cover. Heat will not transfer through multiform topcover, allowing the clinician to hand mould it to the child's foot. Works great off-the-shelf, but customization may be achieved by using a heat gun on the underside of the orthotic.

**SIZING:** Children's shoe size

Size	Children's Shoe Size	Size	Children's Shoe Size
A	4 - 5	E	12 - 13
B	6 - 7	F	1 - 2
C	8 - 9	G	3 - 4
D	10 - 11	H	5 - 6





## P32 Pediatric Hinged Knee Brace

The P32 Pediatric Hinged Knee Brace is designed to provide medial and lateral support to the knee, stabilize the patella, and provide compression and warmth. The hinges may be removed and molded to the leg to provide a custom fit.

**SIZING:**

Size	Knee Centre Circumference
S	12" - 13"
M	13" - 14"
L	14" - 15"



## P34 Pediatric Jumper's Knee Brace

The Pediatric Jumper's Knee Brace is designed to provide compression, warmth, and support to the knee. The inferior and superior straps reduce tension on the patellar tendons, thus reducing the pain associated with repetitive loading of the knee joint.

**SIZING:**

Size	Knee Centre Circumference
S	12" - 13"
M	13" - 14"
L	14" - 15"



## P97 Pediatric Contoured Wrist Brace

The P97 Pediatric Contoured Wrist Support is designed to provide warmth, compression, and stability to the affected wrist. The stays are anatomically contoured to provide a low profile fit that does not impede use of the hand or fingers. Specify left or right.

**SIZING:**

Size	Wrist Circumference
S	3" - 4"
M	4" - 5"
L	5" - 6"

## Custom Order Form - A Knee Supports

### Hinges / Stabilizers

- Polycentric hinges (32 style)
- Multistop hinges (39 style)
- Spiral stays (34 Style)  1  2  3

### Straps

- Elastic circumferential (32 Style -1 superior, 1 inferior)
- Velcro medial to lateral (34 Style - 1 superior, 1 inferior)
- Velcro circumferential (39 style - 2 superior, 2 inferior)

### Buttress Type

- Felt  Removable
- Tube  Donut
- None (patella hole only)  Sewn in superior
- Lateral J  Split J

### Material

- Neoprene
- Coolcel (breathable)
- Airflex (moisture wicking)


### Additional Options

- Pull Tabs
- Popliteal diamond cut out
- Other: \_\_\_\_\_

Notes: \_\_\_\_\_

### Company Information

Company Name : \_\_\_\_\_ Account #: \_\_\_\_\_  
 Patient Name : \_\_\_\_\_  
 Contact Name : \_\_\_\_\_ Phone # \_\_\_\_\_  
 Order Date : \_\_\_\_\_ P.O. # \_\_\_\_\_



### Closure Type

- Regular - closed top and bottom
- Wrap around closure
- Open top closure
- Open bottom closure

### Thickness and Colour

(Neoprene only - colours other than black while supplies last)

- 2mm  Black
- 3mm  Grey
- 6mm  Blue

## Custom Order Form - B Elbow Supports

Please check all boxes that apply

### Hinges

- Free motion hinge     Multistop hinge  
 Bilateral     Medial only     Lateral only

### Material

- Neoprene  
 Coolcel (breathable)  
 Airflex (moisture wicking)

### Thickness and Colour

(Neoprene only - colours other than black while supplies last)

- 2mm     Black  
 3mm     Grey  
 6mm     Blue

### Additional Options

- Hole at elbow  
 Donut support at elbow  
 Tendon pad (universally positionable)  
 Open above elbow closure  
 Open below elbow closure  
 Wrap around closure  
 Other: \_\_\_\_\_

Notes: \_\_\_\_\_

### Company Information

Company Name : \_\_\_\_\_ Account #: \_\_\_\_\_

Patient Name : \_\_\_\_\_

Contact Name : \_\_\_\_\_ Phone # \_\_\_\_\_

Order Date : \_\_\_\_\_ P.O. # \_\_\_\_\_

Please mark all measurements in inches. Please supply two additional measurement half way up and down the brace if total brace length exceeds 10"

Distance from elbow to distal end of brace (normally 5")

Distance from elbow to proximal end of brace (normally 5")

Forearm Circumference

Elbow Centre Circumference

Upper Arm Circumference

Total Length

Left Arm

Right Arm

### Straps

- One adjustable strap above elbow  
 One adjustable strap below elbow  
 Hyperextension cross straps

## Custom Order Form - C Wrist Supports

Please check all boxes that apply

### Stays

- 6" Contoured palmar stay (97A style)
- 6" Contoured dorsal stay (97A style)
- 6" Regular palmar stay (100A style)
- 6" Regular dorsal stay (100A style)
- 8" Long palmar stay (95A style)
- 8" Long dorsal stay (95A style)

### Material

- Neoprene
- Coolcel (breathable)
- Airflex (moisture wicking)

### Thickness and Colour

(Neoprene only - colours other than black while supplies last)

- 2mm       Black
- 3mm       Grey
- 6mm       Blue

### Length

- Standard (6" - 7")
- Long (10" - 11")

Notes: \_\_\_\_\_

### Company Information

Company Name : \_\_\_\_\_ Account #: \_\_\_\_\_

Patient Name : \_\_\_\_\_

Contact Name : \_\_\_\_\_ Phone # \_\_\_\_\_

Order Date : \_\_\_\_\_ P.O. # \_\_\_\_\_

Please mark all measurements in inches. Please supply one additional measurement half way between the wrist and the forearm end of the brace if total brace length exceeds 7"

Palm Circumference

Wrist Circumference

Forearm Circumference

Distance from middle of palm to proximal end of brace

Total Length

Left Arm

Right Arm

### Straps

- 2 adjustable straps      D-ring  Yes  No
- 3 adjustable straps      D-ring  Yes  No
- 4 adjustable straps      D-ring  Yes  No
- Other requests: \_\_\_\_\_

## Custom Order Form - D Shoulder Supports

Please check all boxes that apply

### Material

- Neoprene
- Coolcel (breathable)
- Airflex (moisture wicking)
- Low stretch Neoprene (blue/black only)

### Colour

(Neoprene only - colours other than black while supplies last)

- Black/Red
- Grey/Blue
- Black/Black

### Straps

- Removable webbing strap (UBL Neoprene or Coolcel only)
- Other: \_\_\_\_\_

### Additional Options

- "Sharkskin" non-slip treatment at elbow
- Open at neck closure
- Open at elbow closure
- Wrap around elbow closure
- Torso included (6" unless otherwise stated)
- Other: \_\_\_\_\_

### Company Information


Company Name : \_\_\_\_\_ Account #: \_\_\_\_\_

Patient Name : \_\_\_\_\_

Contact Name : \_\_\_\_\_ Phone # \_\_\_\_\_

Order Date : \_\_\_\_\_ P.O. # \_\_\_\_\_

Right shoulder  Left shoulder



### Measurements

**A** Length from shoulder to elbow

**B** Circumference around shoulder and under armpit

**C** Circumference around chest

**D** Circumference at apex of biceps

**E** Circumference at elbow

Phone: 604-520-3414  
 Toll Free Phone: 1-800-663-1254  
 Fax: 604-520-1193  
 Email: orders@orthoactive.com

# PRODUCT ORDER FORM



Account # \_\_\_\_\_ Phone # \_\_\_\_\_  
 Company Name \_\_\_\_\_ Purchaser's Name \_\_\_\_\_  
 PO # \_\_\_\_\_ Shipping Method  Ground  Air

NOTES: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

ITEM #	LT/RT		SIZE							ITEM #	LT/RT		SIZE						
	M/F		XS	S	M	L	XL	XXL	R		U	M/F		XS	S	M	L	XL	XXL
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			
_____																			

### FOOT PRODUCTS

	A	B	C	D	E	F	G	H	I	1	2	3	4	5	6	7
_____																
_____																
_____																

### MISCELLANEOUS

ITEM #	DESCRIPTION & SIZE
_____	_____
_____	_____
_____	_____